



Compact High Intensity Focused Ultrasound Device ABS Portable HIFU Machine For Home Use

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: Guangzhou, China
- Brand Name: LZR
- Model Number: HIFU
- Minimum Order Quantity: 1 Sets
- Price: Discussible
- Packaging Details: Regular Packaging
- Delivery Time: 5-8 working Days
- Payment Terms: L/C, T/T,
- Supply Ability: 10000 pieces/month




Product Specification

- Name: Portable HIFU Machine For Home Use
- Spacing: 1.0-10mm(0.5mm/step)
- Suitable Skin Types: All
- Specification: Packing , 64 X 81 X 131 Cm
- Material: ABS
- Pack Size: 60*44*46cm
- Lamp Longevity: 300,000 Shots (IPL)
- Scope: 5.0mm-25mm
- Working Mode: Simultaneous / Independent Control
- Highlight: **ABS Portable HIFU Machine , Compact High Intensity Focused Ultrasound Device , Portable HIFU Machine For Home Use**



More Images



Product Description

Compact ABS Home Use HIFU Beauty Device

Product Description:

The HIFU Beauty Machine is made from high-quality ABS material, which ensures that it is durable and long-lasting. This Beauty Therapy Machine also comes with a Lamp Longevity of 300,000 Shots (IPL), which means that it is designed to last for a very long time.

The HIFU Beauty Machine has a vertical body type, which makes it easy to use and store. It also comes with a specification of Packing, 64 X 81 X 131 Cm, which means that it is compact enough to be used in a variety of different settings.

This Beauty Device is suitable for all skin types, which means that anyone can use it to achieve younger and healthier-looking skin.

Whether you have oily, dry, or combination skin, the HIFU Beauty Machine can help improve the look and feel of your skin.

Overall, the HIFU Beauty Machine is an excellent choice for anyone looking for an effective and non-invasive way to improve their skin health. With its high-quality construction, Lamp Longevity of 300,000 Shots (IPL), and suitability for all skin types, this Beauty Therapy Machine is sure to provide you with the results you're looking for.

Fea



tures:

Product Name: HIFU Beauty Machine

Lamp Longevity: 300,000 Shots (IPL)

Length: 5.0 - 25mm (1.0mm/Step)

Pack Size: 60*44*46cm

Specification: Packing, 64 X 81 X 131 Cm
Working Mode: Simultaneous / Independent Control

Our HIFU Beauty Machine is a powerful ultrasound anti-aging device that can help you achieve a more youthful appearance. With its hydrafacial beauty machine capabilities, it can deeply cleanse and exfoliate your skin, leaving it looking fresh and radiant. The Hifu ultrasound facelift machine can also help lift and tighten your skin, reducing the appearance of wrinkles and fine lines. Its long lamp longevity of 300,000 shots (IPL) ensures that it will last for many years to come.

Technical Parameters:

Lamp Longevity	300,000 Shots (IPL)
Size	Portable
Energy Output	0.1J-2.5J (0.1J/step)
Length	5.0 - 25mm (1.0mm / Step)
Energy Type	HIFU (High Intensity Focused Ultrasound)
Body Type	Vertical
Supplier Info	Manufacturer
Spacing	1.0-10mm(0.5mm/step)
Working Mode	Simultaneous / Independent Control
Pack Size	60*44*46cm

This HIFU Face Lift Device is a type of Hifu Ultrasound Facelift Machine that utilizes High Intensity Focused Ultrasound (HIFU) technology to provide non-invasive skin tightening and lifting. It also functions as a Hydrafacial Beauty Machine, promoting healthier-looking skin.





Applications:

The LZR HIFU Beauty Machine is a high-quality Beauty Therapy Machine that is made in Guangzhou, China. The product is designed to provide a non-invasive, safe, and effective way to rejuvenate and refresh the skin. The device can be used in a variety of settings and scenarios, including the following:

Skin Tightening: The LZR HIFU Anti-wrinkle Device is particularly effective in tightening the skin on the face, neck, and other parts of the body. The device is capable of delivering energy output levels ranging from 0.1J-2.5J, making it highly versatile and adaptable to different skin types and conditions. The device can be used to lift and tighten sagging skin, reduce the appearance of fine lines and wrinkles, and improve overall skin tone and texture.

Beauty Salons: The LZR HIFU Beauty Machine is ideal for use in beauty salons and spas. The vertical body type and length of 5.0 - 25mm (1.0mm / Step) make it easy to use and maneuver, even in tight spaces. The device can be used to offer clients a range of skin rejuvenation and tightening treatments, helping them to look and feel their best.

Personal Use: The LZR HIFU Beauty Machine is also suitable for personal use at home. The device is easy to operate and comes with clear instructions, making it a convenient and effective way to maintain healthy, youthful-looking skin. The device can be used to treat a variety of skin concerns, including fine lines, wrinkles, and sagging skin.

Medical Spas: The LZR HIFU Beauty Machine is also well-suited for use in medical spas and clinics. The device has been certified as ISO compliant, ensuring that it meets high standards of quality and safety. The device can be used to offer patients a range of non-invasive skin rejuvenation and tightening treatments, helping them to achieve their desired results without the need for surgery or downtime.

The LZR HIFU Beauty Machine is available for purchase with a minimum order quantity of 1 set. The price of the machine is discussible, and it comes with regular packaging. The delivery time is 5-8 working days, and payment terms include L/C and T/T. The device has a supply ability of 10,000 pieces per month, making it readily available for use in a variety of settings and scenarios.

Overall, the LZR HIFU Beauty Machine is a versatile and effective tool for anyone looking to improve their skin's appearance and health. Whether you are a beauty professional or someone looking to maintain healthy, youthful-looking skin at home, this device is an excellent choice.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou