



LB-100 Home Use Beauty Machine 755nm 808nm 1064nm For Beauty Salon

Our Product Introduction

For more products please visit us on hifubeauty-machine.com


Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: T9
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Home Use Beauty Machine
- Usage: Beauty Salon, Spa, Home
- Output Power: 40W
- Wavelength Range: 430/480/530/560/640/690 ~ 1200nm
- Model: LB-100
- Style: Portable
- Pulse Width: 0.02~0.99s
- Energy Density: 1-10J/cm²
- Application: For Commercial
- Highlight: Home Use Beauty Machine 755nm,
Home Use Beauty Machine 808nm,
Beauty Salon Laser Machine 1064nm



More Images





Product Description

Professional Multi Wavelength LB-100 Home Use Beauty Machine For Salon

Product Description:

Equipped with a pulse width of 0.02~0.99s, the LB-100 Laser Beauty Machine can effectively target a wide range of skin concerns. Whether you're looking to reduce wrinkles, remove unwanted hair, or rejuvenate your skin, this machine has got you covered.

The LB-100 Laser Beauty Machine features a wavelength range of 430/480/530/560/640/690 ~ 1200nm, allowing for precise targeting of different skin types and conditions. This makes it an ideal choice for those seeking a non-invasive alternative to treatments such as Hydrafacial Beauty Machine or Co2 Laser Beauty Machine.

With its three laser bar numbers of 3*755nm, 4*808nm, and 3*1064nm, the LB-100 Laser Beauty Machine can provide effective laser hair removal for a variety of skin types and hair colors. This makes it an excellent home laser hair removal machine for those looking for a safe and convenient way to remove unwanted hair.

Overall, the LB-100 Laser Beauty Machine is a top-of-the-line beauty device that offers a range of benefits to users. Whether you're a beauty professional or an individual seeking to enhance your skincare routine, this machine is definitely worth considering.

Features:

Product Name: Laser Beauty Machine

Usage: Beauty Salon, Spa, Home

Energy Density: 1-10J/cm²

Treatment Heads: 5 Pieces

Laser Bar Number: 3*755nm; 4*808nm; 3*1064nm

Pulse Width: 0.02~0.99s

This Beauty Machine is perfect for use in Beauty Salons, Spas, or even at home. With an Energy Density range of 1-10J/cm² and 5 Treatment Heads, this Home Laser Hair Removal Machine is highly effective. The Laser Bar Number is 3*755nm, 4*808nm, and 3*1064nm, ensuring optimal results for all skin types. The Pulse Width ranges from 0.02~0.99s, making this Home Beauty Machine a versatile and valuable addition to any beauty routine.

Technical Parameters:

Frequency:	1-6Hz
Energy Density:	1-10J/cm ²
Wavelength Range:	430/480/530/560/640/690 ~ 1200nm
Usage:	Beauty Salon, Spa, Home
Treatment Heads:	5 Pieces
Style:	Portable
Screen:	8-inch Touchscreen
Model:	LB-100
Pulse Width:	0.02~0.99s
Treatment Probes:	5 Pieces

Applications:

Equipped with 5 treatment heads, the Laser Beauty Machine can perform a range of functions, from rejuvenation to hair removal. The hydrafacial beauty machine head provides a deep cleansing effect, removing dirt and impurities from the skin. The IPL laser beauty machine head is perfect for reducing the appearance of blemishes, dark spots, and fine lines. Lastly, the home laser hair removal machine head allows users to achieve smooth and hair-free skin in the comfort of their own homes.

Featuring an energy density of 1-10J/cm² and a wavelength range of 430/480/530/560/640/690 ~ 1200nm, the Laser Beauty Machine can effectively target different skin types and colors. The 5 treatment probes can also be interchanged according to the user's needs.

For those looking to improve their skin's texture and appearance, the Laser Beauty Machine is the perfect accompaniment for a regular beauty routine. It can be used in salons or at home, providing a convenient and cost-effective option for beauty enthusiasts. Whether you're looking to reduce the appearance of wrinkles, eliminate unwanted hair, or achieve a clearer complexion, the Laser Beauty Machine has got you covered.

Overall, the Laser Beauty Machine is a must-have for anyone looking to level up their beauty game. With its versatility, convenience, and effectiveness, it's no wonder why it's become a popular choice among beauty professionals and enthusiasts alike.







Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou