



50W Home Use Radio Frequency Beauty Machine 4MHZ With 5 Treatment Tips

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: D101
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



Product Specification

- Name: Radio Frequency Beauty Machine
- Shockwave Frequency: 1-18Hz
- Voltage: 110V/220V10% 60HZ 10A
- Delay Interval: 0-3s
- Rf Applicator: Bipolar
- Suitable For: Beauty Clinic/centre/salon/home Use
- Target Area: Body/Face/Legs/Arms
- Frequency: 4MHZ
- Treatment Tips: 5pcs(F-S/F-M/F-L/I-M/I-L)
- Highlight: 50W Radio Frequency Beauty Machine, Home Use Radio Frequency Beauty Machine, Radio Frequency Beauty Machine 4MHZ



More Images





Product Description

50W Home Use Radio Frequency Beauty Machine With 5 Treatment Tips

Product Description:

Whether you're looking to remove unwanted hair, reduce the appearance of fine lines and wrinkles, or rejuvenate your skin, this facial beauty device has got you covered. It features a delay interval of 0-3 seconds, giving you the flexibility to customize your treatment according to your specific needs.

Designed to operate on both 110V and 220V 10% 60HZ 10A voltage, this beauty machine is suitable for use in a wide range of settings, including spas, beauty salons, and even at home. Its IC equipment technology ensures that it performs reliably and consistently, so you can enjoy the benefits of your beauty treatments with complete peace of mind.

With its advanced laser, E-light, and IPL technology, this beauty device is capable of delivering exceptional results. It can help to minimize the appearance of acne scars, age spots, and other blemishes, leaving your skin looking smoother, brighter, and more youthful. So if you're looking for a powerful and effective beauty device that can help you achieve your desired results, look no further than the RF Beauty Equipment. Order yours today and start enjoying the benefits of beautiful, radiant skin!

Features:

Product Name: RF Beauty Equipment

Rf Power: 50W

Frequency: 4MHZ

Shockwave Frequency: 1-18Hz

Action: IC Equipment

Suitable For: Beauty Clinic/centre/salon/home Use

Beauty Care Instruments

Home Use Beauty Machine

Home Beauty Equipment

Technical Parameters:

Product Parameter

Frequency 50-60HZ

Probe diameter 37mm

Rated voltage 110-220V

Power 72W

Cool temperature -15°

Hot temperature 40°

Technology Cool + hot + EMS + LED

Warranty 1 year





Applications:

This IPL laser beauty machine offers a range of features that make it a must-have for anyone looking to enhance their skin's appearance. With a delay interval of 0-3s, this machine is designed to provide controlled and precise treatment, ensuring maximum safety and efficacy.

Operating at a frequency of 4MHZ, this IPL beauty equipment uses monopolar RF polar technology to deliver deep and effective treatments that target specific areas of the body, face, legs, and arms. Whether you're looking to tighten loose skin, reduce cellulite, or smooth out fine lines and wrinkles, this machine can help you achieve your beauty goals.

Our RF Beauty Equipment is suitable for use in a variety of settings, from professional beauty clinics and centers to salons and even at home. So whether you're a beauty professional looking to expand your range of services, or a beauty enthusiast looking to enhance your own beauty routine, this machine is the perfect choice for you.

So why wait? Invest in our RF Beauty Equipment today and experience the amazing benefits of this powerful and effective IPL beauty equipment for yourself!

Customization:

Our RF Beauty Equipment is an IC equipment that can be customized to meet your specific needs. We offer product customization services that allow you to choose from a wide range of options, including multi-language support, a big color touch screen, and a 4MHZ frequency. This allows you to personalize your device to provide the best user experience for your customers.

Our product is optimized for beauty salons, spas, and clinics, and is suitable for a wide range of skincare treatments. With our E Light Beauty Machine and IPL Laser Beauty Machine, you can offer your customers a variety of skincare treatments that are tailored to their individual needs.

FAQ:

Q: What is RF beauty equipment?

A: RF beauty equipment is a non-invasive cosmetic tool that uses radio frequency energy to tighten and tone the skin, reduce the appearance of wrinkles, and stimulate collagen production.

Q: What are the benefits of using RF beauty equipment?

A: The benefits of using RF beauty equipment include smoother, firmer skin, a reduction in the appearance of wrinkles and fine lines, improved skin texture and tone, and increased collagen production for long-lasting results.

Q: Where is the RF beauty equipment manufactured?

A: The RF beauty equipment is manufactured in Guangzhou, China.

Q: Is RF beauty equipment safe?

A: Yes, RF beauty equipment is a safe and non-invasive cosmetic tool. It is approved by the FDA and has been tested for safety and effectiveness.

Q: How often should I use RF beauty equipment?

A: It is recommended to use RF beauty equipment once or twice a week for best results. However, the frequency of use may vary depending on individual needs and preferences.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou