



## LZR 50W RF Face Beauty Equipment Multi Language With Microwave Shockwave

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: D101
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



#### Product Specification

- Name: Face Beauty Equipment
- Suitable For: Beauty Clinic/centre/salon/home Use
- Language: Multi-language
- Delay Interval: 0-3s
- Action: IC Equipment
- Keywords: Laser E-light Ipl Machine
- Shockwave Frequency: 1-18Hz
- Target Area: Body/Face/Legs/Arms
- Rf Power: 50W
- Highlight: Face Beauty Equipment Multi Language, 50W RF Face Beauty Equipment



#### More Images





## Product Description

### Multi Language Home Beauty Machine Face Beauty Equipment With Microwave Shockwave And 50W RF Power

#### Product Description:

Our RF Beauty Equipment is a multi-language device, making it easy to use no matter where you are in the world. Plus, with its bipolar RF applicator, you can rest assured that you're getting top-of-the-line technology that will provide effective and long-lasting results. Whether you're looking for a beauty device for personal use or for your business, our HIFU RF Machine is the perfect choice. So why wait? Invest in the best and start seeing the results you've been dreaming of!

#### Features:

Product Name: RF Beauty Equipment

Language: Multi-language

Rf Power: 50W

Frequency: 4MHZ

Rf Applicator: Bipolar

Target Area: Body/Face/Legs/Arms

Features:

Aesthetic Medicine Equipment

CET RET RF Machine

IPL Laser Beauty Machine

#### Technical Parameters:

Delay Interval:	0-3s
Treatment Tips:	5pcs(F-S/F-M/F-L/I-M/I-L)
Suitable For:	Beauty Clinic/centre/salon/home Use
Frequency:	4MHZ
Rf Power:	50W
Keywords:	Laser E-light Ipl Machine
Shockwave Frequency:	1-18Hz
Rf Polar:	Monopolar
Screen:	Big Color Touch Screen
Voltage:	110V/220V10% 60HZ 10A

This Home Use Beauty Machine is an IPL Laser Beauty Machine that can be used in Beauty Clinics, Centre, Salons, and also at home.

#### Applications:

One of the key features of this IPL laser beauty machine is its monopolar RF polar technology. This technology is designed to penetrate deep into the skin, helping to stimulate collagen production and tighten the skin. With a powerful 50W RF power input, this device can provide powerful and efficient treatments that can help to deliver real results.

Our IC equipment is designed to be easy to use, with a big color touch screen that provides an intuitive interface for controlling the device. Whether you are a professional beauty technician or just someone who wants to improve their own skin, this machine is easy to use and delivers reliable results.

Some of the key product application occasions and scenarios where this device can be used include:

Beauty salons and spas

Home use by individuals who want to improve their skin

Medical clinics that offer skin treatments

Skin care clinics that focus on facial treatments

This RF beauty equipment product is made in China Guangzhou, which is known for its high-quality manufacturing and attention to detail. With its advanced technology and powerful features, it is a must-have for anyone who is serious about improving their skin's appearance.

#### Customization:

Place of Origin: China Guangzhou  
Voltage: 110V/220V10% 60HZ 10A  
Rf Polar: Monopolar  
Action: IC Equipment  
Delay Interval: 0-3s  
Language: Multi-language

Our Beauty Apparatus is perfect for home use and can be customized to meet your needs. Contact us today to learn more about our product customization services.

## FAQ:

**Q: What is RF beauty equipment?**

A: RF beauty equipment is a non-invasive skin rejuvenation treatment that uses radio frequency energy to improve the appearance of fine lines, wrinkles, and sagging skin.

**Q: Where is the RF beauty equipment made?**

A: The RF beauty equipment is made in Guangzhou, China.

**Q: Is RF beauty equipment safe?**

A: Yes, RF beauty equipment is a safe and non-invasive treatment that has been used for many years in the cosmetic industry.

**Q: How many sessions are needed to see results from RF beauty equipment?**

A: The number of sessions needed to see results from RF beauty equipment varies depending on the individual's skin condition. However, most people see results after 4-6 sessions.

**Q: How long does each RF beauty equipment session take?**

A: Each RF beauty equipment session typically takes between 30-60 minutes depending on the area being treated.





**Guangzhou Lizhirun Technology Development Co., Ltd**

☎ 18688905627

✉ cenglinlin25@gmail.com

🌐 hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou