



## Handheld Introduction Device Microcrystalline Pen Stainless Steel Material

Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: WZ-3
- Minimum Order Quantity: 1
- Price: USD/60-100
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



### Product Specification

- Name: Microcrystalline Pen
- Package Contents: Microneedle Therapy Machine, Power Adapter, User Manual
- Treatment Area: Face, Body
- Speed Levels: 5
- Function: Skin Rejuvenation, Scar Removal, Wrinkle Removal
- Power Supply: AC 100-240V, 50/60Hz
- Usage: Home, Salon
- Needle Length: 0.25mm-2.5mm
- Power Consumption: 10W
- Highlight: Handheld Microcrystalline Pen, Introduction Device Microcrystalline Pen



### More Images





## Product Parameters



### Specification

Brand	Hydra pen	Model	H5
EMS intensity level	4 levels	EMS frequency	0.25 Hz
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Handle Material	ABS
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Battery capacity	700mAh	G.W.	19.4*14.6*6.9 cm/1 pc/ 690 g
Available Catridges	H12, H36, H24, Nano-HS, Nano-HR	1 carton	41*32*38 cm/20 pcs/? kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② Inductive base ③ Hydra pen H5 ④ Cable ⑤ Adapter ⑥ Manual ⑦ Retaining ring* 2 pcs ⑧ 12 pins * 2 pcs

## 4-pole release micro current

Using 4-pole micro current to release molecules that penetrate the bottom layer.



## 6-speed deep penetration

1-6 speed intelligent speed regulation solves various skin problems.



## LED display screen intelligent pause

LED digital display can be used for charging and plugging, making operation more convenient.



## Adjustable needle length 0~1.0 mm

Deep adjustment, The functions of regulation, repair, and hydration are fully equipped.





# EMS | infusion carrier essence device





**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou