



## 140A Derma Stamp 0.25mm-2.5mm Microneedle RF Machine Overheating Protection

### Our Product Introduction

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: WZ-3
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



#### Product Specification

- Name: Microneedle RF Machine
- Needle Length: 0.25mm-2.5mm
- Usage: Home, Salon
- Speed Levels: 5
- Power Supply: AC 100-240V, 50/60Hz
- Model: MTM-01
- Treatment Area: Face, Body
- Color: White
- Needle Number: 12
- Highlight: 0.25mm Microneedle RF Machine,  
2.5mm Microneedle RF Machine



#### More Images





## Product Description

### 140A Derma Stamp Microneedling Face Roller Massager Adjustable Titanium Needle Length Skincare Beard Growth Scalp Hair

#### Product Description:

One of the standout features of the MTM-01 is its safety features. The machine comes equipped with an auto-stop function that prevents the needles from over-penetrating the skin, which can cause damage and discomfort. Additionally, the machine has built-in overheating protection to ensure that it doesn't get too hot during use, which could also damage the skin.

The MTM-01 is made from high-quality stainless steel, which is both durable and easy to clean. The needles themselves have an adjustable length, ranging from 0.25mm to 2.5mm, which allows you to customize your treatment to suit your specific needs and skin type.

In addition to its microneedle therapy capabilities, the MTM-01 can also be used as a butt vacuum therapy machine. Simply attach the vacuum cup to the wand and use it to help tone and firm the skin on your buttocks.

The MTM-01 is powered by an AC 100-240V, 50/60Hz power supply, which means that it can be used in a variety of different countries. Whether you're at home or traveling abroad, you can take your microneedle therapy machine with you and continue your skin health routine no matter where you are.

Overall, the MTM-01 Microneedle Therapy Machine is a versatile and effective tool for anyone looking to improve their skin health and appearance. Its gold RF microneedle technology, adjustable needle length, and safety features make it a standout product in the world of skin care.

#### Features:

Product Name: Microneedle Therapy Machine

Material: Stainless Steel

Package Contents:

Microneedle Therapy Machine

Power Adapter

User Manual

Color: White

Usage: Home, Salon

Safety Features:

Auto-stop

Overheating Protection

#### Technical Parameters:

Speed Levels	5
Material	Stainless Steel
Usage	Home, Salon
Needle Number	12
Color	White
Needle Length	0.25mm-2.5mm
Safety Features	Auto-stop, Overheating Protection
Function	Skin Rejuvenation, Scar Removal, Wrinkle Removal
Package Contents	Microneedle Therapy Machine, Power Adapter, User Manual
Treatment Area	Face, Body

#### Applications:

This Microneedle Therapy Machine is perfect for individuals who want to enhance their skin's texture and appearance. The device is versatile and can be used on any part of the body, including the face, body, and even the scalp. The machine features 12 needles that are

designed to penetrate the skin gently, allowing for maximum absorption of skin care products. The needle length can be adjusted from 0.25mm to 2.5mm, depending on the user's preference and the area being treated. The Microneedle Therapy Machine is an excellent investment for anyone who wants to achieve smoother, more radiant skin. The device is easy to use and can be used at home or in a professional setting. It is perfect for people who are looking for an alternative to invasive skin treatments such as laser therapy or surgery. The Microneedle Therapy Machine is an excellent option for individuals who want to improve the appearance of fine lines, wrinkles, and scars. In addition to its skin care benefits, the Microneedle Therapy Machine can also be used as an Electric Wand Massager. The device features powerful electromagnetic therapy that can help relieve muscle tension and promote relaxation. It is perfect for individuals who suffer from chronic pain or have a busy lifestyle that can cause stress and tension. The Microneedle Therapy Machine can also be used as a Butt Vacuum Therapy Machine. The device features a unique vacuum technology that can help lift and tone the buttocks. It is perfect for individuals who want to achieve a more toned and shapely posterior without the need for invasive surgery or expensive treatments. The Microneedle Therapy Machine comes in a sleek white color that is both stylish and modern. The device is compact and easy to store, making it perfect for individuals who have limited space. It is also durable and long-lasting, ensuring that users get the most out of their investment.





# EMS annular micro current+carrier essence

Continuous output of circular multi-polar current, reaching deep into the skin, resulting in dual efficacy of EMS+essence

## EMS Level 4 Microcurrent

Stimulate cells and accelerate collagen synthesis.



## 4-pole micro current enhances firm skin texture

EMS circular quadrupole energy field. With the aqueous conductivity of essence, micro current and essence can penetrate into the deep skin and achieve great treatment effect.



## Solve kinds of skin problems

- ✓ Skin discoloration
- ✓ Fine lines, wrinkles
- ✓ Large areas of acne
- ✓ Pore size, acne scars
- ✓ Stretch mark on belly
- ✓ Stretch mark on thighs
- ✓ Helps improve skin texture



Red bloodshot




Reduce acne marks



Reduce scars

## Product Parameters




### ● Specification

Brand	Hydra pen	Model	H5
EMS intensity level	4 levels	EMS frequency	0.25 Hz
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Handle Material	ABS
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Battery capacity	700mAh	G.W.	19.4*14.6*6.9 cm/1 pc/ 690 g
Available Catridges	H12, H36, H24, Nano-HS, Nano-HR	1 carton	41*32*38 cm/20 pcs/? kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box   ② Inductive base   ③ Hydra pen H5 ④ Cable   ⑤ Adapter   ⑥ Manual ⑦ Retaining ring* 2 pcs   ⑧ 12 pins * 2 pcs



## Adjustable liquid output

please adjust the discharge of liquid base on the consistency of liquid.



1-3 step fluid adjustment



## Serum Container:



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou