




## Home Spa Facial Microneedle Therapy Machine With Auto Stop Overheating Protection

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)


#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: WZ-3
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



#### Product Specification

- Name: Microneedle Therapy Machine
- Color: White
- Treatment Area: Face, Body
- Power Supply: AC 100-240V, 50/60Hz
- Package Contents: Microneedle Therapy Machine, Power Adapter, User Manual
- Safety Features: Auto-stop, Overheating Protection
- Needle Length: 0.25mm-2.5mm
- Speed Levels: 5
- Model: MTM-01
- Highlight: Home Microneedle Therapy Machine, Facial Microneedle Therapy Machine, Spa Microneedle Therapy Machine



## Product Description

### Home Spa Facial Microneedle Therapy Machine With Auto Stop Overheating Protection

#### Product Description:

The MTM-01 offers unparalleled safety features to ensure a safe experience for the user. The device is equipped with an auto-stop function that stops the device if it detects any abnormal behavior. Additionally, the MTM-01 has overheating protection, which prevents the device from overheating and causing any damage to the skin. These safety features make the MTM-01 a reliable and safe option for anyone looking to improve their skin's appearance.

The MTM-01 is capable of treating various areas of the body, including the face and body. This device is perfect for those looking to improve their skin's texture, reduce wrinkles, and promote collagen growth. The MTM-01 uses RF microneedling technology to penetrate the skin's surface and promote collagen growth in the treated area. This technology is highly effective and provides long-lasting results. The MTM-01 also includes a cold therapy machine that can be used to reduce inflammation and redness after treatment. This feature is particularly useful for those with sensitive skin who may experience discomfort after a treatment session. The cold therapy machine is easy to use and provides instant relief to the treated area.

The MTM-01 also includes a butt vacuum therapy machine that can be used to enhance the appearance of the buttocks. This feature is perfect for those looking to tone and shape their buttocks without having to undergo surgery. The butt vacuum therapy machine is highly effective and can provide noticeable results in just a few sessions.

In conclusion, the Microneedle Therapy Machine (MTM-01) is an excellent device for anyone looking to improve their skin's appearance. This device is safe, effective, and easy to use, making it suitable for both home and salon settings. With its RF microneedling technology, cold therapy machine, and butt vacuum therapy machine, the MTM-01 is a versatile device that can provide a wide range of skin treatments.

#### Features:

Product Name: Microneedle Therapy Machine

Safety Features: Auto-stop, Overheating Protection

Speed Levels: 5

Usage: Home, Salon

Treatment Area: Face, Body

Power Consumption: 10W

This Beauty Massager Machine is perfect for RF microneedling and can be used as an Electric Wand Massager.

#### Technical Parameters:

Usage	Home, Salon
Function	Skin Rejuvenation, Scar Removal, Wrinkle Removal
Treatment Area	Face, Body
Model	MTM-01
Safety Features	Auto-stop, Overheating Protection
Speed Levels	5
Needle Length	0.25mm-2.5mm
Power Consumption	10W
Material	Stainless Steel
Power Supply	AC 100-240V, 50/60Hz

#### Applications:

One of the most popular applications for the WZ-3 is as a facial treatment tool. By gently puncturing the skin with its microneedles, the WZ-3 can stimulate collagen production, reduce fine lines and wrinkles, and improve overall skin texture. When used in conjunction with a BIO Microcurrent Machine, the WZ-3 can also help to lift and tone facial muscles, giving you a more youthful and refreshed appearance. The WZ-3 is also highly effective as a body treatment tool. Its 12 needles can be used to target specific areas of the body, such as stretch marks, scars, and cellulite. Additionally, the WZ-3 can be used in combination with a Beauty Massager Machine to provide a relaxing massage while simultaneously stimulating collagen production and improving skin texture.

For those looking to soothe and calm irritated skin, the WZ-3 can be used as a Cold Therapy Machine. Simply place the device in the refrigerator for a few hours before use, and then use the chilled needles to gently massage the skin. This can help to reduce inflammation, soothe redness, and provide an overall cooling and refreshing sensation.

With a power consumption of just 10W, the WZ-3 is an energy-efficient and cost-effective beauty tool that is perfect for both home and salon use. Whether you're looking to improve your skin's health and appearance or provide your clients with a high-quality beauty treatment, the WZ-3 Microneedle Therapy Machine is the perfect choice.

**Customization:**


---

Place of Origin: China Guangzhou  
Model: MTM-01  
Usage: Home, Salon  
Treatment Area: Face, Body  
Package Contents: Microneedle Therapy Machine, Power Adapter, User Manual  
Speed Levels: 5  
Upgrade to include additional features such as Butt Vacuum Therapy Machine, BIO Microcurrent Machine, or Beauty Massager Machine  
Contact us today to customize your Microneedle Therapy Machine and enhance your skincare routine!



## 4-pole release micro current

Using 4-pole micro current to release molecules that penetrate the bottom layer.



## 6-speed deep penetration

1-6 speed intelligent speed regulation solves various skin problems.



## LED display screen intelligent pause

LED digital display can be used for charging and plugging, making operation more convenient.




## Adjustable needle length 0~1.0 mm

Deep adjustment, The functions of regulation, repair, and hydration are fully equipped.



## Product Parameters




### Specification

Brand	Hydra pen	Model	H5
EMS intensity level	4 levels	EMS frequency	0.25 Hz
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Handle Material	ABS
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Battery capacity	700mAh	G.W.	19.4*14.6*6.9 cm/1 pc/ 690 g
Available Catridges	H12, H36, H24, Nano-HS, Nano-HR	1 carton	41*32*38 cm/20 pcs/2 kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box   ② Inductive base   ③ Hydra pen H5 ④ Cable   ⑤ Adapter   ⑥ Manual ⑦ Retaining ring* 2 pcs   ⑧ 12 pins * 2 pcs

### Adjustable liquid output

please adjust the discharge of liquid base on the consistency of liquid.



### 1-3 step fluid adjustment



### Serum Container:





## Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou