




7D HIFU Beauty Machine Liposonic Vaginal Tightening Gold RF 3 In 1 Beauty Machine

Our Product Introduction

for more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: china
- Brand Name: LZR
- Model Number: 33-9
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Flight case
- Delivery Time: 5-10working days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 300 sets/month



Product Specification

- Name: HIFU Beauty Machine
- Product Name Special 2.0 Ultra Slim Cartridge Face Lifting Body Slimming 7D Hifu Machine Technology: Special 2.0 Ultra Slim Cartridge Face Lifting Body Slimming 7D Hifu Machine Technology: High Focused Ultrasound Cartridge Option: 1.5/2.0/3.0/4.5/6.0/9.0/13.0mm Advantage: Micro+macro Double Treatment Shots: 20000 Shots/ Cartridge Frequency:
- Input Voltage: AC100V-240V
- Treatment Areas: Face, Neck, Body
- Cartridge Shot Number: 10000 Shots
- Type: High Intensity Focused Ultrasound
- Hifu Frequency: 10MHz, 4MHz, 7MHz
- V Max: 60,000 Shots/cartridge
- Highlight: **7D HIFU Beauty Machine, Vaginal Tightening HIFU Beauty Machine, Gold RF 3 In 1 Beauty Machine**



More Images



Product Description

V-Max 7D Hifu Beauty Machine Liposonic Vaginal Tightening Skin Therapy Gold RF 3 In 1 Beauty Device

What is 7D HIFU and why all the craze?

* 7D HIFU Ultrasound is back to serve the beauty. This time, the goal is to restore skin tension.

* The Ultraforms III is an optimum, non-invasive ultrasound device that does skin tightening and facial lifts for a youthful complexion as well as for the body for slimmer contours. With precision per shot pulse, HIFU-powered transducers are designed to either remodel collagen to remove facial wrinkles and sagging skin or tighten body tissues to reflect your true potential.

* Until now, only radiofrequency was proposed in this indication. Radiofrequency devices raise the skin temperature to 45 ° C and only on the outer layer of the skin. We obtain only an ephemeral tensor effect of a few hours related to transient edema. Thus, radiofrequency is inefficient for skin resurfacing.

* That's why we use High Intensity Focused Ultrasound (HIFU), which warms tissues at 70 ° in one point like a magnifying glass under the sun.

Face

Comparable to utilizing sunlight and reflecting it through a magnifying glass to deliver a beam, the ultraforms 3 treatment delivers ultrasound energy into the target region of the skin, stimulating and formulating collagen regeneration to give a skin tightening effect. The precision in delivery of a high density of power naturally triggers neo-collagenesis within the body.

Body

Stubborn fat cells are difficult to burn off even through strenuous exercises and extensive weight loss regimens. The Ultraforms 3 treatment provides both non-invasive and innovative solutions without including side effects or damage to skin.

Advantages:

It uses cartridges concentrated on various depths of skin for lifting and tightening the face just as repairing body tissues for contouring purposes. The patient can get customized treatments on face and body with the steady conveyance of ultrasound energy, quick application, and exact accuracy of MMFU-controlled transducers.

Advantages for patients:

- Applicable for all skin types and tones
- 15–20-minute duration with no downtime
- No pre or post-treatment required
- Quick and durable outcomes

Functions:

Full face lift (overall tightening)
V-shape face lift
Double chin reduction
Jowl reduction (tightens and reduces saggy jowls)
Brow lift
Neck lift











Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebianshan Village, Baiyun District,
Guangzhou