




## Beauty Salon Full Face Skin Test Machine With UV / RGB / PL Light Multilanguage Support

### Our Product Introduction

#### Basic Information

- Place of Origin: china
- Brand Name: Bitmoji
- Certification: CE
- Model Number: A6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month




#### Product Specification

- Name: Face Skin Test Machine
- Function: Skin Test Skin Analysis
- Resolution: Industrial Grade 36 Megapixels
- Feature: RGB White Light, Positively Polarized Light, Negatively Polarized Light, Wu's Light, UV Light, Red Light, Brown Light, Mixed Light
- Report Items: Skin Pore, Wrinkle, Roughness,etc
- High Light: Skin Testing Machine, Skin Analysis Equipment
- Highlight: Beauty Salon Face Skin Test Machine, Full Face Skin Test Machine, Salon Face Skin Analysis Machine



#### More Images



## Product Description

### Beauty Salon Full Face Skin Tester Machine With UV / RGB / PL Light Multilanguage Support

Full face Skin Scanner Analyzer Facial Tester Skin Diagnosis for Salon Use

#### Features:

Epidermal skin analysis, as color of depth, property of skins and other skin problems.

Dermal skin analysis, giving advice on potential skin disorders on dermal layer.

Wrinkle analysis: Property, density, depth and cleaning ways.

Speckles analysis: Analyze depth area, sizes, density accurately and how to dispel of speckles

12 - 15 Mega pixels with auto-focus

Auto-analysis and report mechanism

UV Voice system

Adjustable chine-holder

Multilanguage support

Integrated customer management software and crystal report software

Electrical watchdog protection

Computer and deck are optional fittings

The BS-3200 facial skin analysis system utilizes skin care photographic imaging tool that provides clinical measurement of surface and subsurface facial skin conditions, using digital images technologies via RGB visible light, PL polarized light, and UV spectrum. In use by cosmetic surgeons, dermatologists, and other medical and skin care professionals, the BS-3200 system is a self-contained unit that captures multi-spectral photos of the face, examining the complexion for wrinkles, spots, pores, texture, porphyrins (evidence of bacteria in pores), UV spots (characteristic of photodamage, typically from overexposure to sun), pigmented spots and acne. Additionally, it grades an individual's skin features relative to others of the same age and skin type. Application:

1. Facial skin anlysis to find skin disorders
2. Tell skin potential problems
3. Track skinchange trend
4. Sell beauty products and servi





Bitmoji® 摩玑  
领先科技 创造美丽

A6外贸版目前可支持**16**种语言



中文简体/繁体



俄语



韩语



日语



印尼语



英语



德语



法语



西班牙语



意大利语



土耳其语



荷兰语



越南语



泰语



葡萄牙语



波兰语

更多语言  
敬请期待!



Look forward to more languages!

Bitmoji® 摩玑

THE ACTION REPORT SHOWS  
DIFFERENT VISUAL EXPERIENCE

### Skin problem case show

Your skin diagnostician!

- Multi-dimensional photography
- 8 spectral imaging
- multiple skin testing modes
- 15 detection indicators
- multiple reporting modes
- multiple comparative analysis
- Simulated real beauty skin simulation
- data cloud storage
- File management
- skin test results adjustment
- background system management



— A6-AI intelligent skin tester —

Your skin diagnostician! Specializes in skin management



# Bitmoji®

## A6-AI intelligent skin tester



It currently supports **16** languages

**More language versions are coming soon!**



*Look forward to more languages!*

### A6-AI Intelligent Image Instrument Product Specifications

Your skin diagnostic expert! Professional Skin Care Management Solution



Device Name	A6-AI intelligent skin tester
Model	ZMLH-A6
Pixel	Industrial grade 36 Million pixel
Spectrum	RGB white light; Positive polarized light; Negative polarized light; Wood's light; UV light; Red light; Brown light; Mixed light
Shading Method	Semi-open
Product Material	Industrial Grade ABS
Power Supply	AC 100-260V, 50/60Hz
Product Dimensions	Unfolded: 636 x 521 x 489mm; Folded: 361 x 521 x 489mm
Central Processor	Rockchip RK3568 ; Quad-core Cortex-A55
Motherboard	R10-S6810 Pmotherboard, running Android 11
Main Memory	Dual-channel LP DD R4, 4G
Storage	MMC 5.1, 32G
Operating Method	multi-point capacitive touch
Screen Features	Foldable 90 degrees
WIFI	Built-in dual-band WIFI (2.4G, 5G)
HDMI	1 port (optional)
USB	2 ports
Screen Ratio	16:9
Screen Size	15.6 inches
Screen Resolution	1920*1080
Net Weight	9.5 KG
Gross Weight	14.6 KG
Packaging Dimensions	640 x 570 x 605mm
Packaging Material	Corrugated paper + pearl cotton
Other Accessories	Handheld skin analyzer, power cord













**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeaute-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou