




Beauty Salon Cryolipolysis Slimming Machine 1000W For Non Surgical Fat Reduction

Our Product Introduction

for more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: china
- Brand Name: GLM
- Certification: CE
- Model Number: V10
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month




Product Specification

- Name: Cryolipolysis Slimming Machine
- Technology: Vacuum+Cavitation+ RF
- Keywords: Vela Shape Body Slimming Machine
- Application: Beauty Salon, Spa, Skin Management Center,
- Vacuum Pressure: 0-106Kpa
- RF Frequency: 2.44MHZ
- Cavitation Frequency: 40KHZ
- Pulse Width Of Bio: 100us~300us
- Treatment Area: Body Legs Arms Face Eyes
- Power: 1000W
- High Light: Fat Removal Machine, Portable Fat Freezing Machine
- Highlight: **Beauty Salon Cryolipolysis Slimming Machine, Cryolipolysis Slimming Machine 1000W, Fat Reduction Cryo Slim Machine**



More Images



Product Description

State-Of-The-Art Cryolipolysis Slimming Machine For Non-Surgical Fat Reduction

Overview

V100 has combined with the RF, vacuum, outer automatic roller, infrared light, green light and blue light, cavitation, BIO, lipolaser. It is the new updates technology. It also combines with RF and 3 photons (red, green, blue) to generate heat. It can penetrate the skin to a depth of 15mm by mechanical movement, which is far more than the depth of the other devices currently on the market. While the heat is conducted, the mechanical movement is carried out simultaneously and so is the fat. In this way, it makes the skin more flexible. It can improve the figure and the skin surface. Although the weight of the patient does not reduce, the patient will look much slimmer because the whole fat area has been re-assigned.

The Advantages of Instrument V100

RF+ 2* 40KHZ cavitation +vacuum +IR+BIO all in one handle
 The newest hi-tech body contouring systems
 Most effective slimming machine with all popular functions
 Infrared light heats the tissue up to 2 mm depth
 Radio frequency (RF) heats tissue from 5~15mm depth
 Vacuum+/- massage mechanisms enable precise targeting of energy to the tissue.
 Mechanical massage provides lymphatic drainage (dermal and hypodermal)
 Outer automatic roller with discontinuous will make patients more comfortable during treatment
 This machine can treat for whole body: arms, thighs, buttocks, neck, face, eyes, etc.
 Can be used for 3 clients at the same time
 Combined with the velasMOOTH, lipo laser slimming,
 RF sonic cavitation and BIO for body contouring, shaping,
 fat and cellulite reduction, etc.
 Anti-aging, wrinkle removal, skin lifting&tightening, etc.
 3 touch screens with fashion modern design
 6 units of operation handle and 8 pieces of large lipo laser pads.



Technical Parameters

Input voltage
 AC220V, 50HZ/AC110V, 60HZ
 Power consuion
 1000W
 Fuse size
 10A
 Large vacuum roller handle
 Vacuum pressure
 0-100 (KPA)
 Infrared light wavelength
 700nm
 Green light wavelength
 560nm
 Blue light wavelength
 490nm

Roller speed
50r/min
R.F Energy
30W
R.F Frequency
5MHZ
Middle vacuum roller handle
Vacuum pressure
0-100 (KPA)
R.F Energy
20W
Roller speed
30r/min
R.F Frequency
5MHZ
Infrared light wavelength
700nm
Small vacuum facial handle
Vacuum pressure
0-100(KPA)
RF frequency
5MHZ
Infrared light wavelength
700nm
Vacuum eyes pen (same as small vacuum facial handle)
Small 40KHZ cavitation handle
Cavitation frequency
40KHZ
Large 40KHZ cavitation handle
Cavitation frequency
40KHZ
RF frequency
2.44MHZ
BIO frequency
200HZ~500HZ
Pulse width of BIO
100us~300us
Lipo laser pads
Lipo laser wavelength
650nm
Power
150nw



For a variety of scenarios

Multi-functional handle head, used all over the body, suitable for a variety of scenes, beauty and body industry extension of the essential instrument!



Beauty Salon




Plastic surgery
hospital



Beauty studio



• Skin management center



Multi-dimensional effectiveness

Blasting fat increases body metabolism

Face lift

Improves facial relaxation and tightening

Burst fat

High frequency energy waves break down fat cells

Operation handle description

Multiple function handle head, the whole body can be used, suitable for a variety of scenes, beauty industry to broaden the necessary instruments!



Multifunctional integrated head

The probe is a comprehensive operation head, a collection of 40K burst fat, B10 micro current, vacuum negative pressure, laser, RF RF, automatically scan metabolic aging cells, activate new cell division, to achieve the effect of beauty and skin rejuvenation.



Large roller head

The use of negative pressure suction and unique roller design, the skin and the lower skin muscle tissue up, and then through the role of the roller before and after the vacuum suction, the skin and fat tissue suction, knead, roller mechanical aerobic exercise, so that the body more effectively consume the accumulation of fat, to achieve the effect of slimming body.



Medium size roller head

Using negative pressure suction and unique roller design, the body is deeply massaged in waves to promote the metabolism of the skin, so that the skin can produce more collagen and elastic fiber cells, thereby increasing the elasticity of the skin.



Negative pressure head for face

The use of negative pressure on the face sagging and collagen loss and other problems, all concentrated in a small part through the radio frequency and infrared light to tighten and lift the skin, promote the blood circulation of the face, help the face to enhance and tighten, fade fine lines, show the aesthetic outline.



Radiofrequency eye negative pressure pen

Use negative pressure and RF to help blood circulation in the eyes, eliminate eye bags, improve dry fine lines and dark circles, solve the accumulation of excess fat in the eyes, and promote the absorption of nutrients.



Sail for the perfect body

Multi-function handle for both body and face operation

Vertical V10 Shaper

RADIO FREQUENCY BEAUTY APPARATUS

Vera Shaping Series

Multiple functions Handle head body and face can be operated

- ✓ Full body Slimming
- ✓ Increased elasticity
- ✓ Fade orange peel
- ✓ Face lift




All-in-one all-body lift

Multi-function handle, body and face can be operated

PARAMETER

Product parameter



Product name: V10 Upright

Voltage: 110V~220V

Power: 1000W

Fuse: 10A

Big roller speed: 50 revolutions/min

Medium roller speed: 30 revolutions/min

Ultrasonic frequency: 40KHz

Radio frequency: 2.44MHzB10

Frequency: 200Hz~500Hz

B10 Pulse width: 100 microseconds ~300 microseconds

Negative pressure: Okpa~100kpa

Laser wavelength: 650 nm Power: 150 nanowatts





Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou