



## Effective Fat Removal Slimming Cold Laser Therapy Machine 800W 10D Maxlipo Master Laser

### Our Product Introduction

for more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

#### Basic Information

- Place of Origin: china
- Brand Name: LZR
- Certification: CE
- Model Number: A62
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



#### Product Specification

- Product Name: 10d Maxlipo Master Laser Slimming Machine
- Keywords: 10d Maxlipo Master Laser Machine
- Technology: 10D Maxlipo Master Laser
- 10D Maxlipo Master Laser: 1-200mW(Continuous) / 1-300mW(Pulse)
- Wavelength: 532nm(green Light) / 635nm(red Light) For Option
- Number Of Diode: 10pcs
- Power: 800W
- Application: For Commercial, Beauty Salon, Clinic Etc.
- Function: Pain Relief, Knee Arthritis, Therapy Applications, Body Slimming, Skin Tightening, Diminish Inflammation Etc.
- Treatment Area: About 30mm\*800mm For Body
- High Light: 10d Cold Laser Therapy Machine, Maxlipo Master Cold Laser Therapy Machine, Fat Removal Slimming Cold Laser Machine
- Blue Wavelength: 405nm



#### More Images



## Product Description

### Effective Fat Removal Slimming Cold Laser Therapy 10D Maxlipo Master Laser

#### A Brief Introduction of 10D Laser Maxlipo Master Slim System

With the 532nm cold source laser, the 10 adjustable laser lamps target the fat cell beneath the treatment areas, creating a temporary pore within the fat cell, the fatty liquids can drain out of the cell, the volume of the fat cell is reduced. The fatty acids and glycerol will be flushed out naturally through lymphatic metabolism.

#### About 10D Laser Maxlipo Master Slim System Technical Parameters

Product name 10d maxlipo master laser slimming machine

Keywords 10d maxlipo master laser

Fuse 5A

Wavelength 532nm (green light)

Energy Output 1-200mW (Continuous) / 1-300mW (Pulse)

Treatment Area About 30\*800mm

Mode Continuous, Pulse

Intensity Level 1-3 (adjustable)

Rotation Speed Level 1-3 (adjustable)

Screen 12.1 inch Touch Color Screen

Cooling System Air cooled

Power 800W

N.W. / G.W. 39KG/58KG

Package Size 118cm\*58cm\*60cm









Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou