



1200W 600W Powerful Body Hair Removal Machine / Permanent Hair Removal Machine

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: 808 hair removal instrument
- Certification: CE
- Model Number: NGY-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Body Hair Removal Machine
- Machine Power: 1200W
- Power: 600W
- Machine App: Permanent Laser Hair Removal Machines
- Laser Type: Other
- Skin Types: I-VI
- Style: Stationary, Vertical
- Spotsizes: 12mm*12mm
- Package Size: 75*54*54cm3
- Highlight: 600W Body Hair Removal Machine,
Powerful Body Hair Removal Machine,
1200W Permanent Hair Removal Machine



More Images



Product Description

1200W 600W Powerful Body Hair Removal Machine For All Skin Tones

Product Description:

Weighing only 15 Kg, this hair removal laser machine is lightweight and portable. It can be easily moved around your home or salon, making it the perfect solution for anyone who needs a reliable hair removal system on-the-go.

The cooling fan technology of this laser hair removal machine ensures that the skin is protected during the treatment process. With 2pcs of fans, this machine helps to cool the skin while effectively targeting the hair follicles, making the process painless and comfortable for the user.

The package size of this hair removal device is 75*54*54cm3, making it compact and easy to store. The small size also means that it will not take up too much space in your home or salon, making it a great choice for those who have limited space.

This hair removal laser machine uses a unique laser type that effectively targets the hair follicles, leaving the skin smooth and hair-free.

The laser technology is safe and effective, providing long-lasting results that will leave you satisfied with your purchase.

Safety is a top priority with this hair removal device, which is why it comes equipped with skin contact cooling and automatic power cut-off features. The skin contact cooling feature ensures that the skin is protected during the treatment process, while the automatic power cut-off feature ensures that the machine will turn off if it detects any abnormalities in the treatment process.

In summary, the Laser Hair Elimination Machine is the perfect solution for anyone who wants a fast, effective, and pain-free hair removal experience. With its lightweight and portable design, cooling fan technology, and safety features, this device is the ultimate hair removal solution for anyone who wants to achieve smooth and hair-free skin.



Features:

Product Name: Hair Removal Laser Machine
Energy Density: Up To 120 J/cm2
Power: 600W
Display: 15-inch Touch Screen
Cooling Type: TEC+Sapphire Contact Cooling
Spot Size: 12mm*12mm

This Laser Hair Elimination Machine is a professional hair removal device that features a 15-inch touch screen display and can deliver energy density of up to 120 J/cm2. It has a power output of 600W and comes equipped with TEC+Sapphire Contact Cooling to ensure a comfortable experience. The spot size is 12mm*12mm.



Technical Parameters:

Product Attribute	Technical Parameter
Style	Stationary, Vertical
Display	15-inch Touch Screen
Machine Power	1200W
Power	600W
Cooling Type	TEC+Sapphire Contact Cooling
Energy Density	Up To 120 J/cm2
Cooling Fan	2pcs Fans
Laser Type	Other

Skin Types	I-VI
Weight	15 Kg

This Painless Hair Removal System is a Professional Hair Removal Device with advanced technology for effective hair removal.

Applications:

The Hair Removal Laser Machine is a perfect solution for individuals who want to get rid of unwanted hair permanently. It is suitable for both men and women and can be used on various body parts such as the face, arms, legs, underarms, and bikini area. This device is ideal for those who do not want to undergo painful waxing, threading, or shaving sessions.

The Hair Removal Laser Machine comes with a TEC+Sapphire Contact Cooling system that ensures a cooling effect on the skin during the hair removal process, making it comfortable for the user. The spot size of 12mm*12mm ensures that the device can cover a large area and reduces the time taken for the hair removal process.

The Hair Removal Laser Machine is perfect for use in beauty salons, spas, and clinics, providing professional hair removal services to customers. It can also be used at home by individuals who want to enjoy the benefits of permanent hair removal at their convenience.

The package size of 75*54*54cm3 ensures that it is easy to store and transport.

In conclusion, the Hair Removal Laser Machine is a perfect solution for anyone who wants to get rid of unwanted hair permanently. It is a professional hair removal device that utilizes other laser type technology suitable for skin types ranging from I-VI. The TEC+Sapphire Contact Cooling system ensures a comfortable hair removal process. With a spot size of 12mm*12mm, it can cover a large area and reduce the time taken for the hair removal process. It is suitable for use in beauty salons, spas, clinics, and at home.



Customization:

Looking for a painless hair removal system that can help you get rid of unwanted hair? Look no further than our 808 hair removal instrument, model number NGY-1, made in China Guangzhou.

Our laser hair elimination machine features TEC+Sapphire contact cooling technology, which ensures maximum comfort during use. The machine also comes equipped with 2 cooling fans to help regulate temperature and prevent overheating.

Our unwanted hair remover device is perfect for professional use in salons and spas, or for personal use at home. The pulse duration can be adjusted from 10-400ms, providing customizable treatment options for a variety of skin types and hair colors.

Invest in our permanent laser hair removal machine today and enjoy smooth, hair-free skin in no time.





Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou