



## Automated Microneedling Device 5 Levels Microneedling Electric Pen 0.25mm-2.5mm

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: M8
- Minimum Order Quantity: 1
- Price: SUS\$20-40
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



#### Product Specification

- Name: Automated Microneedling Device
- Frequency: 50-60Hz
- Needle Length: 0.25mm-2.5mm
- Needle Material: Medical Grade Stainless Steel
- Input Voltage: 110-240V
- Size: 16.5cm X 4.5cm X 4.5cm
- Usage: Skin Rejuvenation, Acne Scar Treatment, Wrinkle Removal
- Speed Levels: 5 Levels
- Power Source: Electric
- Highlight: Automated Microneedling Device 5 Levels, Microneedling Electric Pen 0.25mm, Microneedling Electric Pen 2.5mm



## Product Description

### Automated Microneedling Device With 5 Adjustable Speed Levels For At Home Skin Renewal

#### Product Description:

Measuring 16.5cm X 4.5cm X 4.5cm, this device is compact and easy to use. With a frequency of 50-60Hz and an output voltage of 5V, it's powerful enough to give you the desired results.

The Electric Micro Needle Device features medical-grade stainless steel needles that are designed to penetrate the skin gently. With nine injections in the hand, it gives you an even and precise application, ensuring that all areas of your skin receive the necessary treatment.

This device also features a crystal cannula 32g, which assists with the application of serums and other skin care treatments. The cannula ensures that the serum is evenly distributed into the skin, providing maximum absorption and better results.

Using the Electric Micro Needle Device is simple and easy. Start by cleaning your face and applying your favorite serum. Turn on the device and adjust the needle depth according to your preference. Move the device gently over your face in a circular motion, focusing on problem areas.

The Electric Micro Needle Device is perfect for anyone looking to improve the texture and appearance of their skin. It's ideal for reducing the appearance of fine lines, wrinkles, acne scars, and hyperpigmentation.

In conclusion, the Electric Micro Needle Device is a must-have for anyone looking to improve their skin's overall health and appearance. With its powerful motor, precise needles, and crystal cannula, this device is the perfect addition to your skincare routine.

#### Features:

Product Name: Electric Microneedle Machine

Needle Material: Medical Grade Stainless Steel

Usage: Skin Rejuvenation, Acne Scar Treatment, Wrinkle Removal

Input Voltage: 110-240V

Color: White

Needle Number: 12 Needles

The Electric Microneedle Machine is an Automated Needle Therapy System designed to provide skin rejuvenation, acne scar treatment, and wrinkle removal. It uses medical grade stainless steel needles and has 12 needles for maximum efficiency. The device operates on an input voltage of 110-240V and is available in white color.

#### Applications:

One of the most popular applications of the Electric Micro Needle Device is in skin rejuvenation treatments. The device can effectively reduce the appearance of fine lines, wrinkles, and scars, leaving the skin looking younger and more radiant. It is also used in the treatment of acne scars and hyperpigmentation, providing a safe and effective solution for these common skin problems.

Another common use for the Electric Micro Needle Device is in hair restoration treatments. The device can be used to stimulate hair growth by increasing blood circulation and promoting the delivery of nutrients to the hair follicles. This can help to improve hair growth and thickness, making it an ideal solution for those suffering from hair loss or thinning hair.

In addition to its beauty applications, the Electric Micro Needle Device can also be used in medical settings. It is commonly used in the treatment of stretch marks, surgical scars, and burns, providing a non-invasive and effective solution for these conditions.

The Electric Micro Needle Device is a lightweight device that weighs only 1.5kg, making it easy to transport and use in a variety of settings. Its sleek white design and model number, M5, make it a stylish addition to any beauty or medical practice. It operates on a wide input voltage range of 110-240V and frequency of 50-60Hz, making it suitable for use in various countries and regions.

In summary, the Electric Micro Needle Device is a versatile and effective tool that can be used in a wide range of applications and scenarios. Its nine injections in the hand and advanced technology make it a popular choice for beauty and medical professionals alike. Whether you are looking to rejuvenate your skin, stimulate hair growth, or treat a medical condition, this device is a safe and effective solution.

## 产品/尺寸

PRODUCT  
SIZE :





**Customization:**



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou