



## Stainless Steel Electric Micro Needle Machine For Facial Rejuvenation

Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: H5
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



### Product Specification

- Name: Electric Micro Needle Machine
- Weight: 1.5kg
- Color: White
- Frequency: 50-60Hz
- Needle Number: 12 Needles
- Output Voltage: 5V
- Needle Material: Medical Grade Stainless Steel
- Input Voltage: 110-240V
- Usage: Skin Rejuvenation, Acne Scar Treatment, Wrinkle Removal
- Highlight: **Stainless Steel Electric Micro Needle Machine, Facial Rejuvenation Electric Micro Needle Machine**



### More Images





## Product Description

### Medical Grade Stainless Steel Microneedle Machine For Facial Rejuvenation

#### Product Description:

The Electric microneedle Machine is powered by electricity and has an output voltage of 5V. This ensures that the device is powerful enough to penetrate the skin's surface to stimulate collagen production, but gentle enough to prevent any damage or discomfort. The device comes with 12 needles that are made of high-quality stainless steel and are arranged in a circular pattern, designed to create micro-channels in the skin when in use. The needles are also replaceable, making it easy to maintain and clean the device after each use. The Electric microneedle Machine also comes with a Crystal Cannula 32g, which is designed to help deliver serums and other topical treatments directly into the skin. This enhances the effectiveness of the treatment and ensures that the skin is receiving the nutrients it needs to rejuvenate and heal.

One of the key benefits of the Electric microneedle Machine is that it can be used for a variety of skin concerns, including acne scar treatment, wrinkle removal, and skin rejuvenation. The device works by creating micro-channels in the skin, which stimulates collagen production and improves the overall texture and appearance of the skin.

In conclusion, the Electric microneedle Machine is an excellent investment for anyone looking to improve their skin's appearance and health. With its Electric Skin Needling Machine technology, Crystal Cannula 32g, and 12 needles, this device is capable of delivering professional-grade results in the comfort of your own home. Whether you're looking to reduce the appearance of acne scars, fine lines, or wrinkles, this device can help you achieve your skincare goals.

#### Features:

Product Name: Electric microneedle Machine

Usage: Skin Rejuvenation, Acne Scar Treatment, Wrinkle Removal

Output Voltage: 5V

Size: 16.5cm X 4.5cm X 4.5cm

Needle Number: 12 Needles

Color: White

Beauty micro crystal pen

De Martha nine needles

Nine injections in the hand

#### Applications:

The Electric microneedle Machine comes in a sleek white color that adds elegance to your skincare routine. It has a needle length range of 0.25mm-2.5mm, which means you can adjust the depth of the needle according to your skin's needs. The device also has five speed levels that allow you to customize the intensity of the treatment.

The De Martha nine needles are perfect for various skin concerns such as fine lines, wrinkles, acne scars, hyperpigmentation, and uneven skin texture. This product is suitable for anyone who wants to improve their skin's overall health and appearance.

The Electric microneedle Machine can be used in various scenarios and occasions. It is perfect for home use, allowing you to achieve professional-level skin treatments in the comfort of your own home. It can also be used in beauty salons, spas, and clinics, providing an effective and safe treatment option for clients.

This product is easy to use and can be incorporated into your skincare routine with ease. It is recommended to use the Electric microneedle Machine once a week, and you can see a noticeable improvement in your skin's texture and appearance after just a few treatments.

In conclusion, the Electric Skin Needling Machine is a versatile and powerful product that can help you achieve your desired skin goals. Its De Martha nine needles, Medical Grade Stainless Steel material, adjustable needle length, and speed levels make it a must-have product in your beauty arsenal. Whether you are a skincare enthusiast or a beauty professional, this product is perfect for you.

## 产品/尺寸

PRODUCT  
SIZE :





## 产品/配件

PRODUCT :  
ACCESSORIES :





### Customization:

Customize your Electric microneedle Machine with our product customization services:

Place of Origin: China Guangzhou

Speed Levels: 5 Levels

Needle Number: 12 Needles

Model: M5

Usage: Skin Rejuvenation, Acne Scar Treatment, Wrinkle Removal

Input Voltage: 110-240V

Our Electric microneedle Machine is a top-of-the-line Dermal Pen System that offers nine injections in the hand for precise and effective treatment. Customize your machine to meet your specific needs and provide your clients with the best possible experience.



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou