




Automated Electric Micro Needle Device For Effortless Skin Rejuvenation

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Hydra.pen
- Certification: CE
- Model Number: H5
- Minimum Order Quantity: 1
- Price: USD/60-100
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Micro Needle Device
- Input Voltage: 110-240V
- Usage: Skin Rejuvenation, Acne Scar Treatment, Wrinkle Removal
- Model: M5
- Size: 16.5cm X 4.5cm X 4.5cm
- Weight: 1.5kg
- Frequency: 50-60Hz
- Needle Material: Medical Grade Stainless Steel
- Output Voltage: 5V
- Highlight: **Automated Micro Needle Device,
Skin Rejuvenation Micro Needle Device,
Electric Micro Needle Device**



Automated Needle Therapy For Effortless Skin Rejuvenation Electric Nano Water Light Introduction Instrument

Product Description:

The Electric Micro Needle Device has a compact size of 16.5cm X 4.5cm X 4.5cm, making it easy to use and maneuver. It is also lightweight, making it an ideal choice for those who want to achieve a more youthful and radiant complexion without the hassle of heavy equipment. The sleek white color of the device adds a touch of elegance to your beauty routine, making it a great addition to any bathroom counter.

This device operates at a frequency of 50-60Hz, making it highly effective at delivering the microcrystal treatment to your skin. The needle material used in this device is medical-grade stainless steel, ensuring that it is safe and hygienic for use on your skin. With its powerful and precise performance, you can trust the Electric Micro Needle Device to deliver the results you desire.

With the Electric Micro Needle Device, you can achieve a number of beauty goals, including reducing the appearance of fine lines and wrinkles, improving skin elasticity, and reducing the appearance of scars and stretch marks. This device is easy to use, and with regular use, you can see visible improvements in your skin's appearance.

In conclusion, the Electric Micro Needle Device is an innovative and effective tool for anyone who wants to improve the appearance of their skin. With its compact size, sleek design, and powerful performance, this device is a must-have for anyone who wants to achieve a more youthful and radiant complexion. So why wait? Order your Electric Micro Needle Device today and experience the benefits for yourself!

Features:

Product Name: Electric microneedle Machine

Output Voltage: 5V

Power Source: Electric

Color: White

Needle Length: 0.25mm-2.5mm

Speed Levels: 5 Levels

Crystal Cannula 32g

Beauty needle

Automated Needle Therapy System

Applications:

The needle material used in this product is medical grade stainless steel, which ensures that the needles are safe and hygienic to use. The device has a needle length range of 0.25mm-2.5mm, which makes it suitable for a range of different treatments and skin types. The device is also lightweight, weighing only 1.5kg, which makes it easy to use for extended periods of time.

The electric microneedle machine is perfect for a range of different beauty treatments. One of the most popular treatments is the Nine Injections in the Hand, which involves using the device to inject hyaluronic acid into the skin on the back of the hand. This treatment helps to reduce the appearance of wrinkles and fine lines, and can also improve the overall texture and appearance of the skin.

Another popular treatment using this device is the Beauty Micro Crystal Pen. This treatment involves using the device to create micro-channels in the skin, which allows for the deeper penetration of serums and other skincare products. This treatment helps to improve the overall appearance of the skin, and can also help to reduce the appearance of acne scars, stretch marks, and other skin imperfections. The electric microneedle machine is a versatile and effective tool that can be used in a range of different scenarios. Whether you are a professional esthetician looking to enhance your services, or an individual looking to take your skincare routine to the next level, this Micro Needle Pen Device is the perfect choice for you.



4-pole release micro current

Using 4-pole micro current to release molecules that penetrate the bottom layer.



6-speed deep penetration

1-6 speed intelligent speed regulation solves various skin problems.



LED display screen intelligent pause

LED digital display can be used for charging and plugging, making operation more convenient.



Adjustable needle length 0~1.0 mm

Deep adjustment,.The functions of regulation, repair, and hydration are fully equipped.



Serum Container:



Adjustable liquid output


please adjust the discharge of liquid base on the consistency of liquid.



1-3 step fluid adjustment



Product Parameters



● Specification


Brand	Hydra pen	Model	H5
EMS intensity level	4 levles	EMS frequency	0.25 Hz
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Handle Material	ABS
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Battery capacity	700mAh	G.W.	19.4*14.6*6.9 cm/1 pc/ 690 g
Available Catridges	H12, H36, H24, Nano-HS, Nano-HR	1 carton	41*32*38 cm/20 pcs/? kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② Inductive base ③ Hydra pen H5 ④ Cable ⑤ Adapter ⑥ Manual ⑦ Retaining ring* 2 pcs ⑧ 12 pins * 2 pcs

EMS annular micro current+carrier essence

Continuous output of circular multi-polar current, reaching deep into the skin, resulting in dual efficacy of EMS+essence

EMS Level 4 Microcurrent

Stimulate cells and accelerate collagen synthesis.



4-pole micro current enhances firm skin texture

EMS circular quadrupole energy field. With the aqueous conductivity of essence, micro current and essence can penetrate into the deep skin and achieve great treatment effect.

Solve kinds of skin problems

- ✓ Skin discoloration
- ✓ Fine lines, wrinkles
- ✓ Large areas of acne
- ✓ Pore size, acne scars
- ✓ Stretch mark on belly
- ✓ Stretch mark on thighs
- ✓ Helps improve skin texture



Red bloodshot



Reduce acne marks



Reduce scars

Triple effect to restore skin vitality

EMS Microcurrent+Carrier Liquid Apply to Facial Skin



EMS microcurrent

promotes skin collagen synthesis



Carrier essence

360° coverage penetration



Intelligent repair

Activate skin repair mechanism



Customization:

Our Electric Microneedle Machine is a state-of-the-art beauty needle system that is designed to provide skin rejuvenation, acne scar treatment, and wrinkle removal. Made of medical-grade stainless steel, this Automated Microneedle System is adjustable from 0.25mm to 2.5mm in needle length and is powered by electricity with an input voltage of 110-240V. Our Dermal Pen System is made in China Guangzhou and is customizable to meet your unique needs.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou