



## 7200K Skin Analysis Camera With RGB / UV Lights CCD1/1.7 Sensor Standard Light

Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)


### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Bitmoji
- Certification: ISO CE
- Model Number: A6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



### Product Specification

- Name: Skin Analysis Camera
- PDF Print: YES
- Operating System: RGB +UV
- USB: 2.0 Port
- Computer Requirement: Windows XP/Vista/7/8/10
- Packing Size: L68xW49xH51cm
- Camera Sensor Model: CCD1/1.7
- Color Temperature: 7200 K
- Power Supply: AC 100-240V, 50/60Hz
- Highlight: 7200K Skin Analysis Camera,  
UV Lights Skin Analysis Camera,  
RGB Lights Skin Analyzer Camera



### More Images





## Product Description

### Skin Analysis Camera With RGB/UV Lights CCD1/1.7 Sensor Standard Light Cosmetic Analysis Device

#### Product Description:

With the ability to detect 12 different skin problems, this skin analyzer is a must-have tool for anyone looking to improve their skin care routine. The USB 2.0 port allows for easy connectivity to your computer, making it simple to transfer and store your skin analysis results. With a magnifier rate of 2X, this skin analyzer allows you to see your skin in greater detail than ever before. It also operates using RGB+UV, providing even more detailed information about your skin's health.

For computer requirements, this skin analyzer is compatible with Windows XP/Vista/7/8/10 operating systems.

Whether you're a professional esthetician or simply looking to improve your at-home skin care routine, this skin analyzer is the perfect tool for you. So why wait? Get your hands on our Handheld Skin Analyzer today and start taking control of your skin's health!

#### Features:

Product Name: Skin Analysis Machine

USB: 2.0 Port

Weight: 2.5kg

Color Temperature: 7200 K

Magnifier Rate: 2X

PDF Print: Yes

Features:

Digital Skin Camera Microscope

Face Recognition Body Temperature Scanner

VISIA

#### Technical Parameters:

Technical Parameter	Description
Product Name	DPL Skin Rejuvenation Machine, Acne Removal Machine, Handheld Skin Analyzer
Camera Sensor Model	CCD1/1.7
Detect Problem	12 Skin Problem Detected
Magnifier	Yes
Color Temperature	7200 K
Power Supply	AC 100-240V, 50/60Hz
Working Environment	10 – 35 Celsius Degrees
Operating System	RGB +UV
Weight	2.5kg
Magnifier Rate	2 X
Computer Requirement	Windows XP/Vista/7/8/10

#### Applications:

One of the main application occasions for the Skin Analysis Machine is in beauty salons and spas. Estheticians and beauty therapists can use this device to perform a professional skin analysis for their clients. The Skin Analysis Machine can provide a detailed report on a client's skin type, moisture levels, pigmentation, and other important factors. Based on this information, the esthetician can recommend personalized skin care treatments and products that will be most effective for the client.

The Skin Analysis Machine is also ideal for use in dermatology clinics. Dermatologists can use this device to assess the overall health of a patient's skin and identify any potential issues or conditions. The Skin Analysis Machine can provide detailed images of the skin, including UV damage, age spots, and other blemishes. This information can be used to develop a treatment plan that is tailored to the patient's specific needs.

The DPL Skin Rejuvenation Machine is another application scenario for the Skin Analysis Machine. This device uses LED light therapy to stimulate collagen production and promote skin rejuvenation. The Skin Analysis Machine can be used to assess the effectiveness of the DPL Skin Rejuvenation Machine treatments by providing before-and-after images and data on the patient's skin condition.

Finally, the Skin Analysis Machine can be used by individuals who are interested in monitoring their own skin health and improving their skin care routine. This device is easy to use and can be connected to a Windows XP/Vista/7/8/10 computer. The Skin Analysis Machine can help individuals identify areas of their skin that need extra attention, such as dry patches or pigmentation. By using the Skin Analysis

Machine regularly, individuals can track their progress and make informed decisions about their skin care routine. The Skin Analysis Machine is made in China Guangzhou and has a Color Temperature of 7200 K. Its Portable design makes it easy to use in a variety of settings. The Packing Size measures L68xW49xH51cm, and it requires an AC 100-240V, 50/60Hz Power Supply.



Chinese simplified/traditional



Russian



Korean



Japanese



Indonesian



English



German



French



Spanish



Italian



Turkish



Dutch



Vietnamese



Thai



Portuguese



Polish

It currently supports **16** languages

More language versions are coming soon!

Look forward to more languages!

### A6-AI Intelligent Image Instrument

#### Product Specifications

Your skin diagnostic expert! Professional Skin Care Management Solution



Device Name	A6-AI intelligent skin tester
Model	ZMLH-A6
Pixel	Industrial grade 36 Million pixel
Spectrum	RGB white light; Positive polarized light; Negative polarized light; Wood's light; UV light; Red light; Brown light; Mixed light
Shading Method	Semi-open
Product Material	Industrial Grade ABS
Power Supply	AC 100-260V, 50/60Hz
Product Dimensions	Unfolded: 636 x 521 x 489mm; Folded: 361 x 521 x 489mm
Central Processor	Rockchip RK3568 ; Quad-core Cortex-A55
Motherboard	R10-S6810 Pmotherboard, running Android 11
Main Memory	Dual-channel LP DD R4, 4G
Storage	MMC 5.1, 32G
Operating Method	multi-point capacitive touch
Screen Features	Foldable 90 degrees
WIFI	Built-in dual-band WIFI (2.4G, 5G)
HDMI	1 port (optional)
USB	2 ports
Screen Ratio	16:9
Screen Size	15.6 inches
Screen Resolution	1920*1080
Net Weight	9.5 KG
Gross Weight	14.6 KG
Packaging Dimensions	640 x 570 x 605mm
Packaging Material	Corrugated paper + pearl cotton
Other Accessories	Handheld skin analyzer, power cord



### Customization:

Customize Your Face Moisturizer Machine with Our Skin Analysis Machine!

Our Facial Skin Analyzer, originating from China Guangzhou, comes with a 2X magnifier rate and a built-in magnifier, allowing you to detect 12 different skin problems with ease. Additionally, our Skin Analysis Machine comes equipped with a 2.0 USB port and the ability to print PDFs, making it a versatile tool for your skincare needs.

With these product customization services, you can create personalized and effective Customized Skin Care Products tailored to your unique needs.



**Guangzhou Lizhirun Technology Development Co., Ltd**

☎ 18688905627

✉ cenglinlin25@gmail.com

🌐 hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebain Village, Baiyun District,  
Guangzhou