



RGB UV Skin Analysis Machine 2X Magnification Camera Face Skin Analyzer For Windows

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: ISO CE
- Model Number: A6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: RGB UV Skin Analysis Machine
- Usb: 2.0 Port
- Camera Sensor Model: CCD1/1.7
- Working Environment: 10 – 35 Celsius Degrees
- Type: Portable
- Color Temperature: 7200 K
- Computer Requirement: Windows XP/Vista/7/8/10
- Detect Problem: 12 Skin Problem Detected
- Packing Size: L68xW49xH51cm
- Highlight: RGB UV Skin Analysis Machine,
ISO UV Skin Analysis Machine,
Magnification Camera Face Skin Analyzer



More Images





Product Description

RGB UV 2X Magnification Camera Skin Analyzer For Windows

Product Description:

One of the best features of the Skin Analysis Machine is its ability to print out PDF reports. This means that you can keep a record of your skin analysis results, which can be incredibly helpful when it comes to tracking your progress over time. The Skin Analysis Machine is also equipped with a color temperature of 7200 K, which ensures that you get accurate color representation in your skin analysis.

The Skin Analysis Machine comes with a magnifier that has a rate of 2X. This allows you to zoom in on specific areas of your skin, making it easier to identify any issues that may be present. The magnifier is a great feature for anyone who wants to take a closer look at their skin and identify any potential problems.

The Skin Analysis Machine is perfect for anyone who wants to take control of their skincare routine. By providing you with detailed information about your skin, you can make informed decisions about the products you use. With this information, you can choose customized skin care products that are tailored to your specific needs, ensuring that you get the best possible results.

Overall, the Skin Analysis Machine is one of the best facial skin care devices on the market. Its accuracy and precision make it a valuable tool for anyone who wants to improve their skin health. Whether you're a professional esthetician or someone who wants to take better care of their skin at home, the Skin Analysis Machine is an excellent investment that will help you achieve your skincare goals.

Features:

Product Name: Skin Analysis Machine

... ..

Weight: 2.5kg

... ..

Magnifier: Yes

... ..

Color Temperature: 7200 K

... ..

Operating System: RGB +UV

... ..

Computer Requirement: Windows XP/Vista/7/8/10

... ..

Use: Skin Analyzer

... ..

Additional Function: Acne Removal Machine

... ..

Technical Parameters:

| Product Attribute | Value |
|----------------------|-------------------------|
| Product Name | Body Analyzer Machine |
| Type | Portable |
| Camera Sensor Model | CCD1/1.7 |
| Color Temperature | 7200 K |
| Megapixel | 2000W |
| Usb Port | 2.0 |
| Magnifier | Yes |
| Working Environment | 10 – 35 Celsius Degrees |
| Computer Requirement | Windows XP/Vista/7/8/10 |
| Weight | 2.5kg |
| Pdf Print | YES |

Applications:

The Skin Analysis Machine is perfect for a variety of occasions and scenarios. It can be used in beauty salons, spas, and dermatology clinics to help professionals diagnose and treat skin conditions. It can also be used by individuals at home to monitor the health of their skin and track their progress over time.

One of the key features of the Skin Analysis Machine is its facial recognition body temperature scanner. This allows it to detect changes in body temperature that may be indicative of a skin condition or underlying health issue. With this information, professionals can provide personalized recommendations to their clients and individuals can take proactive steps to improve their skin health. The Skin Analysis Machine also boasts a high-quality camera sensor model, CCD1/1.7, with a megapixel of 2000W. This provides clear and detailed images of the skin, allowing professionals to see even the smallest details and imperfections. Additionally, the magnifier rate of 2 X helps to further enhance the accuracy of the analysis. Overall, the Skin Analysis Machine is an essential tool for anyone interested in maintaining healthy and beautiful skin. Its versatility and accuracy make it an ideal choice for a wide range of applications, from professional settings to personal use.







A6-AI Intelligent Image Instrument

Product Specifications

Your skin diagnostic expert! Professional Skin Care Management Solution



| Device Name | A6-AI intelligent skin tester |
|----------------------|--|
| Model | ZMLH-A6 |
| Pixel | Industrial grade 36 Million pixel |
| Spectrum | RGB white light; Positive polarized light; Negative polarized light; Wood's light; UV light; Red light; Brown light; Mixed light |
| Shading Method | Semi-open |
| Product Material | Industrial Grade ABS |
| Power Supply | AC 100-260V, 50/60Hz |
| Product Dimensions | Unfolded: 636 x 521 x 489mm; Folded: 361 x 521 x 489mm |
| Central Processor | Rockchip RK3568 ; Quad-core Cortex-A55 |
| Motherboard | R10-S6810 Pmotherboard, running Android 11 |
| Main Memory | Dual-channel LP DD R4, 4G |
| Storage | MMC 5.1, 32G |
| Operating Method | multi-point capacitive touch |
| Screen Features | Foldable 90 degrees |
| WIFI | Built-in dual-band WIFI (2.4G, 5G) |
| HDMI | 1 port (optional) |
| USB | 2 ports |
| Screen Ratio | 16:9 |
| Screen Size | 15.6 inches |
| Screen Resolution | 1920*1080 |
| Net Weight | 9.5 KG |
| Gross Weight | 14.6 KG |
| Packaging Dimensions | 640 x 570 x 605mm |
| Packaging Material | Corrugated paper + pearl cotton |
| Other Accessories | Handheld skin analyzer, power cord |

Customization:

Customize your Facial Skin Analyzer with our product customization services. Our Skin Analysis Machine originates from Guangzhou, China and boasts a Megapixel of 2000W. It's compatible with Windows XP/Vista/7/8/10 operating systems and features RGB+UV technology for accurate skin analysis. The machine comes equipped with a magnifier for a closer look and is packaged in a convenient size of L68xW49xH51cm.

With our customization services, you can add the latest Acne Removal Machine technology to your 3D Skin Analyzer. Give your customers the best skin analysis experience with the ability to remove acne using the same machine. Contact us today to learn more about our customization options.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou