



LZR Portable Diode Laser Hair Removal Machine With Adjustable Pulse Duration

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: NGY-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Portable Diode Laser Hair Removal Machine
- Machine App: Permanent Laser Hair Removal Machines
- Coolingtype: TEC+Sapphire Contact Cooling
- Style: Stationary, Vertical
- Display: 15-inch Touch Screen
- Machine Power: 1200W
- Energy Density: Up To 120 J/cm2
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Power: 600W
- Highlight: Portable Diode Laser Hair Removal Machine, CE Portable Diode Laser Machine



More Images





Product Description

Professional Diode Laser Hair Removal Machine With Adjustable Pulse Duration And TEC Sapphire Cooling

Product Description:

The Hair Removal Laser Machine is a state-of-the-art device that uses laser technology to remove hair. This device emits a high-intensity beam of light that penetrates the skin and targets the hair follicles. The laser energy heats up the hair follicles, damaging them and preventing them from growing back. This process is known as Laser Hair Elimination.

The Hair Removal Laser Machine has a pulse duration of 10-400 Ms, which is adjustable depending on the skin type and the area being treated. The device has a power of 1200W, which ensures that it can effectively remove hair from all parts of the body.

The Hair Removal Laser Machine has a cooling system that is designed to protect the skin from the heat generated by the laser. The device has two cooling fans that help to dissipate the heat, and it uses TEC+Sapphire Contact Cooling to ensure that the skin is not damaged during the hair removal process.

The Hair Removal Laser Machine is designed to be stationary and vertical, which makes it easy to use and maneuver. The device has a user-friendly interface that allows users to adjust the settings according to their needs. The device is also equipped with a safety feature that ensures that the laser does not emit unless it is in contact with the skin.

The Hair Removal Laser Machine is an ideal Unwanted Hair Remover Device that is suitable for use by both men and women. It is a safe and effective way to eliminate unwanted hair from different parts of the body, including the face, arms, legs, underarms, and bikini area.

The device is also suitable for all skin types, including sensitive skin.

In conclusion, the Hair Removal Laser Machine is a Laser Hair Elimination Machine that is designed to provide a safe, fast, and effective solution for hair removal. Its advanced cooling system and adjustable pulse duration make it a reliable choice for permanent hair removal. With its high power and user-friendly interface, the Hair Removal Laser Machine is sure to provide users with a hassle-free and comfortable hair removal experience.

Features:

Product Name: Hair Removal Laser Machine


Machine App: Permanent Laser Hair Removal Machines

Display: 15-inch Touch Screen







Style: Stationary, Vertical

Cooling Fan: 2pcs Fans

Safety Features: Skin Contact Cooling, Automatic Power Cut-off





		
Ice Hair Removal	Selectable Spot Size	Diode Laser
		
Fast Light	Adjustable Frequency	Adjustable Energy

Technical Parameters:

Energy Density	Up To 120 J/cm2
Machine Power	1200W
Cooling Type	TEC+Sapphire Contact Cooling
Pulse Duration	10-400 Ms
Skin Types	I-VI
Machine App	Permanent Laser Hair Removal Machines
Safety Features	Skin Contact Cooling, Automatic Power Cut-off
Style	Stationary, Vertical
Laser Type	Other
Power	600W

Applications:

This Hair Removal Laser Machine is perfect for both personal and professional use, as it is easy to operate and can be used for a wide range of hair removal applications. It is ideal for use in beauty salons, spas, and clinics that offer hair removal services. The machine comes with a 15-inch touch screen display, which makes it easy for users to navigate the settings and adjust the parameters to suit their needs. The stationary, vertical style of the machine allows for easy use and accessibility.

The Hair Removal Laser Machine is a powerful system, with a power output of 600W. This makes it an efficient and fast machine that can remove hair in a short period of time, making it ideal for busy salons and clinics. It is also a cost-effective solution for those who want to remove hair at home, as it requires minimal maintenance and has a long lifespan.

Some of the occasions and scenarios where the Hair Removal Laser Machine can be used include:

Beauty salons – the machine can be used to provide hair removal services to customers who want a painless and effective solution for unwanted hair.

Spas – the machine can be used as part of a hair removal treatment package for customers who want a relaxing and rejuvenating experience.

Clinics – the machine can be used to provide hair removal treatments for patients who want a safe and effective solution for unwanted hair.

Personal use – the machine can be used at home by individuals who want to remove unwanted hair without the cost or hassle of going to a salon or clinic.

The Hair Removal Laser Machine (Model Number: NGY-1) is a versatile and effective solution for painless hair removal. Its safety features, touch screen display, and powerful performance make it an ideal choice for both personal and professional use. Its flexibility and efficiency make it an essential addition to any beauty salon, spa, or clinic that offers hair removal services.







Support and Services:

The Hair Removal Laser Machine is equipped with advanced technology to provide safe and effective hair removal treatment. Our team of technical support experts is available to assist you with any questions or concerns you may have regarding the use of the machine. In addition to our technical support, we also offer a range of services to help you get the most out of your Hair Removal Laser Machine. These services include:

Training sessions to ensure proper usage and safety

Maintenance and repair services to keep your machine in top condition

Access to replacement parts and accessories

Consultation services to help you develop effective treatment plans for your clients

Our goal is to provide you with the highest level of support and services to ensure your satisfaction with the Hair Removal Laser Machine.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou