



Permanent Laser Hair Removal And Tattoo Removal Machine With 2 Cooling Fans

Our Product Introduction

For more products please visit us on hifubeauty-machine.com


Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: NGY-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month




Product Specification

- Name: Laser Hair Removal And Tattoo Removal Machine
- Machine App: Permanent Laser Hair Removal Machines
- Laser Type: Other
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Power: 600W
- Skin Types: I-VI
- Pulse Duration: 10-400 Ms
- Coolingtype: TEC+Sapphire Contact Cooling
- Display: 15-inch Touch Screen
- Highlight: **Laser Hair Removal And Tattoo Removal Machine**



More Images



Product Description

600W Advanced Smooth Skin Laser Machine For Permanent Hair Reduction Tattoo Hair Removal Machine

Product Description:

The Hair Removal Laser Machine comes with 2 cooling fans that help to cool the skin during the treatment process. This ensures that the skin does not get damaged due to the heat generated by the laser. The cooling fans are designed to provide maximum comfort to the user during the hair removal process.

One of the best things about this unwanted hair remover device is that it is suitable for all skin types from I-VI. This means that people with all skin types can use this device to remove unwanted hair. The laser hair treatment device uses TEC+Sapphire Contact Cooling technology to cool the skin during the treatment process. This ensures that the skin does not get damaged due to the heat generated by the laser.

The Hair Removal Laser Machine uses an advanced laser type that is different from traditional laser hair removal devices. This laser type is designed to remove unwanted hair quickly and effectively without causing any damage to the skin. The energy density of this device is up to 120 J/cm², which means that it can remove hair from the root effectively.

Say goodbye to unwanted hair and hello to smooth, hair-free skin with the Hair Removal Laser Machine. This unwanted hair remover device is easy to use and can be used from the comfort of your own home. It is a safe and effective way to get rid of unwanted hair permanently.





Features:

Product Name: Hair Removal Laser Machine
Skin Types: I-VI
Safety Features: Skin Contact Cooling, Automatic Power Cut-off
Laser Type: Other
Machine App: Permanent Laser Hair Removal Machines
Style: Stationary, Vertical
Keywords: Unwanted Hair Remover Device, Laser Hair Treatment Device, Laser Hair Elimination Machine

Technical Parameters:

Safety Features	Skin Contact Cooling, Automatic Power Cut-off
Cooling Type	TEC+Sapphire Contact Cooling
Laser Type	Other
Energy Density	Up To 120 J/cm2
Cooling Fan	2pcs Fans
Skin Types	I-VI
Machine App	Permanent Laser Hair Removal Machines
Pulse Duration	10-400 Ms
Style	Stationary, Vertical
Machine Power	1200W

This Hair Removal Laser Machine is a painless hair removal system, unwanted hair remover device, and permanent hair removal device.



Applications:

The NGY-1 Hair Removal Laser Machine is a professional and effective permanent laser hair removal machine designed for use in a variety of settings. This painless hair removal system is a perfect choice for beauty salons, spas, and clinics as well as for personal use at home. The machine is manufactured in China Guangzhou, ensuring high quality and durability.

One of the main advantages of the NGY-1 Hair Removal Laser Machine is its power. With a 600W laser, it is capable of delivering quick and efficient treatment to a wide range of skin types and hair colors. The machine also features two cooling fans to ensure comfortable treatment and prevent overheating. The TEC and Sapphire Contact Cooling technology also guarantees a safe and painless hair removal experience.

The machine's versatility and effectiveness make it perfect for a variety of hair removal scenarios. It can be used to remove hair from the legs, arms, underarms, bikini line, and even the face. With its permanent hair removal capabilities, users can expect long-lasting results with minimal maintenance required.

Beauty salons and spas can use the NGY-1 Hair Removal Laser Machine to offer their clients a painless and effective hair removal experience. The machine's quick and efficient treatment also makes it ideal for busy salons with a high volume of clients. Additionally, personal use at home is also a great option for those looking for a convenient hair removal solution.

In conclusion, the NGY-1 Hair Removal Laser Machine is a powerful, efficient, and versatile hair removal system that can be used in a variety of settings. Its painless hair removal system and permanent hair removal capabilities make it an excellent choice for personal use at home or for beauty salons and spas looking to provide high-quality hair removal services to their clients.

Support and Services:

The Hair Removal Laser Machine is a high-quality product designed to provide effective and long-lasting hair removal results. Our technical support and services team is dedicated to ensuring that you receive the best possible experience with your machine.

We offer a range of technical support services to help you get the most out of your Hair Removal Laser Machine. Our team is available to answer any questions you may have about the product, including troubleshooting, maintenance, and repair services.

In addition to our technical support services, we also offer a range of training and educational resources to help you use your Hair Removal Laser Machine effectively and safely. Our team of experts can provide you with tips and advice on how to get the best results from your machine, as well as guidance on how to use it safely and responsibly. Whether you need technical support, maintenance services, or training and education, our team is here to help you every step of the way. Contact us today to learn more about our services and how we can help you get the most out of your Hair Removal Laser Machine.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou