



## Stationary Vertical Painless Hair Removal Machine 600W Salon Hair Removal Equipment

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)


#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: NGY-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



#### Product Specification

- Name: Painless Hair Removal Machine
- Pulse Duration: 10-400 Ms
- Cooling Fan: 2pcs Fans
- Laser Type: Other
- Machine Power: 1200W
- Style: Stationary, Vertical
- Skin Types: I-VI
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Power: 600W
- Highlight: Stationary Painless Hair Removal Machine, Vertical Painless Hair Removal Machine, 600W Salon Hair Removal Equipment



## Product Description

### Stationary Vertical 1200W Laser Hair Removal Machine For Skin Types I VI

#### Product Description:

The Hair Removal Laser Machine comes equipped with a 15-inch touch screen display that allows for easy navigation and precise operation. The device features TEC+Sapphire contact cooling technology that ensures a comfortable and safe hair removal experience. The cooling fan system with 2pcs fans further enhances the cooling effect, making the process pain-free and enjoyable.

The Hair Removal Laser Machine features a pulse duration of 10-400 ms, which means that it can be adjusted according to the user's skin type and hair thickness. This feature allows for maximum efficacy and minimal discomfort during the hair removal process. The device is also equipped with a water flow sensor that monitors the water flow rate, ensuring that the machine is always functioning optimally.

The Hair Removal Laser Machine is a highly effective and reliable permanent hair removal device that is suitable for use in professional salons as well as at home. With its advanced technology and painless hair removal system, this machine is the perfect solution for anyone looking for a long-lasting and effective hair removal solution. So, if you're tired of traditional hair removal methods and want a permanent solution, the Hair Removal Laser Machine is the perfect device for you.

#### Features:

Product Name: Hair Removal Laser Machine

Display: 15-inch Touch Screen

Machine Power: 1200W

Cooling Fan: 2pcs Fans

Laser Type: Other

Safety Features: Skin Contact Cooling, Automatic Power Cut-off

This Professional Hair Removal Device is the perfect solution for achieving smooth skin. The Smooth Skin Laser Machine is equipped with a 15-inch Touch Screen and 1200W of power. The Painless Hair Removal System uses 2pcs Fans for cooling and has safety features such as Skin Contact Cooling and Automatic Power Cut-off.

#### Technical Parameters:

Technical Parameter	Description
Style	Stationary, Vertical
Power	600W
Machine Power	1200W
Cooling Type	TEC+Sapphire Contact Cooling
Energy Density	Up To 120 J/cm <sup>2</sup>
Skin Types	I-VI
Safety Features	Skin Contact Cooling, Automatic Power Cut-off
Pulse Duration	10-400 Ms
Machine App	Permanent Laser Hair Removal Machines
Laser Type	Other

This Permanent Hair Removal Device is a painless hair removal system that uses a hair depilator machine to remove hair.



### Applications:

With an energy density of up to 120 J/cm<sup>2</sup>, this Laser Hair Elimination Machine is capable of delivering powerful laser pulses that target hair follicles and prevent hair growth. The cooling system, which combines TEC and Sapphire Contact Cooling, ensures that the skin is protected from damage during the treatment process.

The NGY-1 Laser Hair Elimination Machine is designed for use in a variety of settings, including beauty salons, spas, and clinics. Its stationary, vertical design makes it easy to use and allows for precise targeting of hair follicles on various parts of the body.

This Laser Hair Elimination Machine is suitable for use on a wide range of skin types, from I to VI. Whether you have light or dark skin, this machine can be adjusted to deliver the optimal amount of energy for your skin type.

So, whether you're looking to remove unwanted hair from your legs, arms, underarms, or bikini area, the NGY-1 Laser Hair Elimination Machine is the perfect solution. With its powerful energy output and advanced cooling system, it delivers safe and effective hair removal results every time.

Don't settle for a subpar hair removal solution. Invest in the NGY-1 Laser Hair Elimination Machine today and achieve the smooth, hair-free skin you've always wanted!




# PRODUCT DETAILS

Laser type	Diode laser
Laser wavelength	808nm/755nm/1064nm /940nm
Pulse frequency	1~10Hz
Spot size	12*12mm/12*35mm
Pulse Width	1-400 milliseconds (adjustable)
skin type	I-VI skin types
Laser module power	600W/800W/1000W11200W/1600W
Voltage 电压	110V/220V+20V 50/60HZ
Cooling system	Air + water + compressor cooling
Energy	1-100j

# THE PRINCIPLE

## —— 原则

The Diode laser is operating at a wavelength of Alex755nm, 808nm and 1064nm, 3 different wave-lengths come out at the same time to work in the different depth of hair to work full range permanent hair removal result. Alex755nm delivering powerful energy is absorbed by melanin chromophore, making it ideal for skin type 1, 2 and fine thin hair. The longer wavelength 808nm works deeper hair follicle with less absorption of melanin, which is more safety for dark skin hair removal. 1064nm works as infrared red with high water absorption, it is specialized for dark skin hair removal including tanned skin.



# MULTIPLE OPTIONS





# 755NM

Suitable for blonde and fine hair particularly effective

# 808NM

Golden wavelength suitable for all hair types

# 940NM

More suitable for light skin and hair

# 1064NM

For dark, tanned skin particularly effective

## Support and Services:

The Hair Removal Laser Machine is a high-tech product that provides effective and long-lasting hair removal results. Our team of technical support and service professionals are available to assist you with any questions or concerns you may have regarding the use and maintenance of your Hair Removal Laser Machine.

We offer a range of technical support services, including:

Phone support

Email support

Online chat support

Remote assistance via screen-sharing

In addition to technical support, we also offer a variety of services to ensure the optimal performance and longevity of your Hair Removal Laser Machine:

Regular maintenance and cleaning

Repairs and replacements of parts

Training on proper use and safety precautions

Our goal is to provide you with the best possible experience using your Hair Removal Laser Machine. Please do not hesitate to contact us if you need any assistance.



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeaauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou