



120J/Cm2 Laser Hair Removal Salon Machine 600W Unwanted Hair Removal Machine

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: NGY-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Laser Hair Removal Salon Machine
- Pulse Duration: 10-400 Ms
- Skin Types: I-VI
- Laser Type: Other
- Style: Stationary, Vertical
- Coolingtype: TEC+Sapphire Contact Cooling
- Power: 600W
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Energy Density: Up To 120 J/cm2
- Highlight: 120J/Cm2 Laser Hair Removal Salon Machine, Laser Hair Removal Salon Machine 600W, 600W Unwanted Hair Removal Machine



Product Description

TEC Sapphire Contact Cooling Laser Hair Removal Machine With Up To 120 J/Cm2 Energy Density

Product Description:

With a 15-inch Touch Screen Display, operating the Hair Removal Laser Machine has never been easier. The user-friendly interface allows for quick and easy adjustments to the intensity level and duration of the treatment. The display also provides real-time feedback on the energy level, ensuring that the treatment is both safe and effective.

The Painless Hair Removal System is equipped with advanced safety features such as Skin Contact Cooling and Automatic Power Cut-off. Skin Contact Cooling ensures that the skin is protected and cooled before, during, and after treatment. This reduces the risk of skin damage and makes the treatment more comfortable for the user. The Automatic Power Cut-off feature ensures that the machine stops emitting laser energy when it is not in contact with the skin, which prevents accidental exposure to the laser.

Using the Painless Hair Removal System is a breeze - simply adjust the intensity level and duration of the treatment, and let the Hair Removal Laser Machine do the rest. The laser targets the hair follicles, destroying them at the root without damaging the surrounding skin. This results in a painless and effective hair removal experience.

With the Painless Hair Removal System, you can say goodbye to painful waxing and tedious shaving. This Smooth Skin Laser Machine is designed to make hair removal quick, easy, and painless. So why wait? Get your Painless Hair Removal System today and experience the joy of smooth, hair-free skin!

Features:

Product Name: Hair Removal Laser Machine

Safety Features: Skin Contact Cooling, Automatic Power Cut-off

Energy Density: Up To 120 J/cm2

Style: Stationary, Vertical

Skin Types: I-VI

Power: 600W

This Hair Removal Laser Machine is the perfect Unwanted Hair Remover Device that will leave you with smooth skin. Its powerful 600W motor ensures efficient hair removal with its Energy Density of Up To 120 J/cm2 while the Safety Features, Skin Contact Cooling, and Automatic Power Cut-off ensure a Painless Hair Removal System. The stationary, vertical style makes it easy to use on all Skin Types: I-VI.

Technical Parameters:

Product Attribute	Value
Product Name	Unwanted Hair Remover Device
Machine App	Permanent Laser Hair Removal Machines
Laser Type	Other
Energy Density	Up To 120 J/cm2
Coolingtype	TEC+Sapphire Contact Cooling
Cooling Fan	2pcs Fans
Safety Features	Skin Contact Cooling, Automatic Power Cut-off
Power	600W
Machine Power	1200W
Skin Types	I-VI
Display	15-inch Touch Screen

three-wavelength

755NM+808NM+1064NM+940NM, Suitable for all hair colors and skin types
适合所有发色和肤色

high-power

Two 1200W/1600W/2000W high-power handle options available for OEM/ODM
两个1200W/1600W/2000W大功率手柄可供OEM/ODM

Spot size

12*35cm and 12*35cm
12 * 35厘米and12 * 35厘米



12*35cm



PRODUCT DETAILS

Laser type	Diode laser
Laser wavelength	808nm/755nm/1064nm /940nm
Pulse frequency	1~10Hz
Spot size	12*12mm/12*35mm
Pulse Width	1-400 milliseconds (adjustable)
skin type	I-VI skin types
Laser module power	600W/800W/1000W11200W/1600W
Voltage 电压	110V/220V+20V 50/60HZ
Cooling system	Air + water + compressor cooling
Energy	1-100j

755NM

Suitable for blonde and fine hair particularly effective

808NM

Golden wavelength suitable for all hair types

940NM

More suitable for light skin and hair

1064NM

For dark, tanned skin particularly effective

Applications:

This permanent hair removal device is perfect for individuals who have unwanted hair growth on their legs, arms, underarms, and other areas of the body. It is also ideal for those who are tired of the regular waxing and shaving routines that only provide a temporary solution. NGY-1 Hair Removal Laser Machine is designed to provide a long-lasting solution to hair removal problems, which is why it is called a permanent hair removal device.

The cooling fan of the NGY-1 Hair Removal Laser Machine is designed with two fans that help to prevent the machine from overheating. The safety features of the machine are also top-notch, as it is designed with skin contact cooling and automatic power cut-off to prevent any form of skin damage.

This permanent hair removal device has an energy density of up to 120 J/cm², which makes it an excellent option for users who have thick hair growth. The laser hair treatment device is also perfect for individuals who are looking to save money in the long run as they won't have to keep paying for waxing or shaving services regularly.

The NGY-1 Hair Removal Laser Machine is perfect for use in various scenarios, including homes, beauty salons, spas, and other beauty establishments. It is an excellent option for individuals who want to look good at all times, without having to go through the stress of regular hair removal routines.

Support and Services:

The Hair Removal Laser Machine is a high-tech beauty device that uses laser technology to remove unwanted hair in a safe and effective manner. We offer comprehensive technical support and services for this product, including:

- Online user manuals and instructional videos to help you get the most out of your device
- Troubleshooting guides to address any issues you may encounter
- Access to our team of experienced technicians who are available to assist with any technical questions or concerns
- Warranty coverage to ensure that your device is protected against any defects or malfunctions
- Optional maintenance and repair services to keep your device in top working condition

With our commitment to customer satisfaction and reliable support, you can trust that your Hair Removal Laser Machine will provide you with a long-lasting and enjoyable experience.





18688905627



cenglinlin25@gmail.com



hifubeaity-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou