



5 In 1 Home Beauty Machine , Beauty Salon / Home Laser Hair Removal Machine

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: Q9
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: 5 In 1 Home Beauty Machine
- Energy Density: 1-10J/cm²
- Frequency: 1-6Hz
- Pulse Width: 0.02~0.99s
- Treatment Heads: 5 Pieces
- Screen: 8-inch Touchscreen
- Output Power: 40W
- Treatment Probes: 5 Pieces
- Usage: Beauty Salon, Spa, Home
- Highlight: 5 In 1 Home Beauty Machine,
5 In 1 Beauty Salon Hair Removal Machine,
CE Home Laser Hair Removal Machine



More Images



Product Description

5-In-1 Home Beauty Machine Beauty Salon Machine For Professional Spa Treatments

Product Description:

This beauty device is perfect for commercial use and can be used to treat a wide range of skin concerns, such as hair removal, skin rejuvenation, pigmentation, and acne. It is built to provide a safe and comfortable treatment experience, making it an ideal choice for both home beauty machine users and professionals.

If you are looking for a reliable and efficient beauty device, the Laser Beauty Machine is the perfect choice. It is designed to cater to all your beauty needs and provide a safe and comfortable treatment experience. Invest in this machine today and achieve your desired beauty goals!

Features:

Product Name: IPL Laser Beauty Machine

Treatment Heads: 5 Pieces

Screen: 8-inch Touchscreen

Energy Density: 1-10J/cm²


Wavelength Range: 430/480/530/560/640/690 ~ 1200nm

Spot Diameter: 6mm-14mm Adjustable

Technical Parameters:

IPL Laser Beauty Machine	
E Light Beauty Machine	
Home Laser Hair Removal Machine	
Model	LB-100
Screen Size	8-inch Touchscreen
Frequency	1-6Hz
Treatment Probes	5 Pieces
Treatment Heads	5 Pieces
Output Power	40W
Energy Density	1-10J/cm ²
Wavelength Range	430/480/530/560/640/690 ~ 1200nm
Application	For Commercial







Applications:

With its pulse width of 0.02~0.99s, output power of 40W, and energy density of 1-10J/cm², this Diode Laser Beauty Machine is a powerful tool that can treat a variety of skin conditions. It comes equipped with 3*755nm, 4*808nm, and 3*1064nm laser bar numbers that can be used for different skin types and conditions.

The Laser Beauty Machine can be used in a variety of occasions and scenarios. Here are some of the most common ones:

Beauty Salons: Beauty salons can use the Laser Beauty Machine to offer a range of services to their clients, including hair removal, skin rejuvenation, and acne treatment. The portable design of the device makes it easy to move around the salon and use it on different clients.

Home Beauty Machine: The Laser Beauty Machine is also suitable for home use. If you're someone who prefers to take care of their skin at home, this device can help you achieve salon-like results. It's easy to use and comes with instructions that can help you get started.

Medical Clinics: Medical clinics can use the Laser Beauty Machine to treat a variety of skin conditions, including rosacea, spider veins, and pigmentation disorders. The device is safe and effective when used by trained professionals.

Spas: Spas can use the Laser Beauty Machine to offer a range of services to their clients, including hair removal, skin tightening, and wrinkle reduction. The device is versatile and can be used on different parts of the body.

The Laser Beauty Machine is a versatile device that can be used in different occasions and scenarios. Whether you're a beauty salon owner, a medical professional, or someone who loves taking care of their skin at home, this Diode Laser Beauty Machine can help you achieve your skin goals.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaity-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou