



2000M Energy Picosecond Laser Machine 1064nm 532nm For Tattoo Removal / Acne Scars

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: P6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Picosecond Laser Machine
- Screen: 12inch
- Aiming Beam: Diode 655 Nm (Red), Adjustable Off/on
- Rf Frequency: 2M
- Wave Length: 1064nm, 532nm
- Touch Screen Size: 12 Inch
- Voltage: AC220v±10% 50Hz, 110v±10% 60Hz
- Laser Pulse: Pico Laser
- Maximum Energy: 2000M
- Highlight: **2000M Energy Picosecond Laser Machine,
Picosecond Laser Machine 1064nm,
Picosecond Tattoo Removal Machine 532nm**



Product Description

2000M Picosecond Laser Machine With Energy Density ≤ 0.6 J/Cm2 And Spot Size 2-10 Mm

Product Description:

The Picosecond Laser Machine is designed to operate on both AC220v $\pm 10\%$ 50Hz and 110v $\pm 10\%$ 60Hz power supply systems, making it highly versatile and suitable for use in different settings. With a maximum energy output of 2000M, this laser machine is capable of delivering high-intensity laser pulses that can effectively treat various skin conditions.

The Picosecond Laser Machine is equipped with a highly advanced ND:YAG Solid Laser light source, which is renowned for its exceptional performance and reliability. This light source is designed to deliver highly precise laser pulses that are capable of penetrating deep into the skin, making it ideal for treating a wide range of skin conditions, including Fractional infrared laser skin tightening.

The Picosecond Laser Machine comes with a range of features that make it easy to use and operate. One of the standout features of this laser machine is its adjustable spot size, which can be set between 2-10mm. This feature allows for highly precise targeting of specific areas of the skin, making it ideal for treating a wide range of skin conditions.

In conclusion, the Picosecond Laser Machine is a highly advanced and versatile Aesthetic Laser Machine that is designed to help individuals achieve their ideal look. With its advanced Pico Laser technology, high-energy output, and adjustable spot size, this laser machine is capable of delivering highly precise laser pulses that can effectively treat a wide range of skin conditions, including Fractional infrared laser skin tightening.

Features:

Product Name: Picosecond Laser Machine

Spot Size: 2-10 Mm

Rf Frequency: 2M

Other Name: Picosecond Laser

Screen: 12inch

Voltage: AC220v $\pm 10\%$ 50Hz, 110v $\pm 10\%$ 60Hz

Features:

Fractional infrared laser skin tightening

Laser beauty machine

Laser skin rejuvenation machine



Technical Parameters:

Technical Parameters	Description
Voltage	AC220v±10% 50Hz, 110v±10% 60Hz
Aiming Beam	Diode 655 Nm (Red), Adjustable Off/on
Other Name	Picosercond Laser
Touch Screen Size	12 Inch
Wave Length	1064nm, 532nm
Screen	12inch
Energy Density	≤ 0.6 J/cm ²
Laser Pulse	Pico Laser
Maximum Energy	2000M
Spot Size	2-10 Mm

This table shows the technical parameters of the Picosercond Laser machine, also known as the Fractional Laser Machine or Aesthetic Laser Machine. It has a touch screen size of 12 inches and a screen of 12 inches. The laser has a wavelength of 1064nm and 532nm and an energy density of ≤ 0.6 J/cm². The laser pulse is Pico Laser and the maximum energy is 2000M. The aiming beam is a diode 655 Nm (Red) and adjustable Off/on. The spot size is 2-10 Mm, and it is suitable for Fractional infrared laser skin tightening.



Applications:

The LZR P6 Picosecond Laser Machine is perfect for use in various beauty centers, spas, and dermatology clinics. Its advanced technology allows for various skin treatments like skin rejuvenation, tattoo removal, pigmentation removal, acne scars removal, and many more. It is also known as the Superpicosecond freckle machine, which is an excellent solution for people who want to get rid of freckles and other skin blemishes.

One of the best features of the LZR P6 Picosecond Laser Machine is its 12-inch screen that displays all the necessary information clearly. This makes it easier for doctors and aestheticians to perform the required treatments with precision and accuracy. The laser pulse technology used in this machine is Pico Laser, which ensures that the laser penetrates the skin deeper, making the treatment more effective.

The LZR P6 Picosecond Laser Machine is suitable for various skin types, and it comes with adjustable voltage, which makes it easier to use in different countries. It runs on $AC220V \pm 10\%$ 50Hz, $110V \pm 10\%$ 60Hz voltage, making it highly efficient and effective.

Overall, the LZR P6 Picosecond Laser Machine is an excellent investment for anyone in the beauty industry. It is a highly advanced 810 skin rejuvenation machine that offers various skin treatments, making it a versatile machine that can handle multiple skin issues. Its efficiency and effectiveness make it a great tool for any doctor or aestheticians looking to provide high-quality skin treatments to their clients.

Customization:

Customize your LZR P6 Picosecond Laser Machine to fit your specific needs with our product customization services! Our laser skin rejuvenation machine is designed with the latest technology and features a ND:YAG solid laser light source with a maximum energy of 2000M. Our P6 model is manufactured in China Guangzhou and has an RF frequency of 2M with an energy density of $\leq 0.6 \text{ J/cm}^2$. With voltage options of $AC220V \pm 10\%$ 50Hz and $110V \pm 10\%$ 60Hz, our product can be used worldwide.

Our team offers product customization services to ensure that your LZR P6 Picosecond Laser Machine is tailored to your unique needs.

Whether you need additional energy output or a different voltage option, we can modify your machine to meet your requirements. With our services, you can be confident that you are receiving a high-quality aesthetic laser machine that is customized to your specifications. Invest in the best with our LZR P6 Picosecond Laser Machine. Our product is a top-of-the-line IPL Laser Beauty Machine that is designed to provide exceptional results. With our customization services, we can help you achieve your desired outcomes with a product that is tailored to your needs.



Support and Services:

The Picosecond Laser Machine is a cutting-edge technology that offers precise and effective laser treatments for various cosmetic and medical purposes. Our product technical support and services are designed to ensure that our customers receive the highest level of satisfaction with their purchase.

Our technical support team is available to answer any questions or concerns you may have about the operation or maintenance of your Picosecond Laser Machine. We offer comprehensive training to ensure that you and your staff are fully equipped to operate the machine safely and effectively.

In addition to technical support, we also provide a range of services to help you get the most out of your Picosecond Laser Machine. Our service offerings include:

Installation and setup

Regular maintenance and cleaning

Calibration and testing

Repair and replacement of parts

Upgrades and software updates

We understand the importance of keeping your Picosecond Laser Machine in top working condition, which is why we offer prompt and reliable service to minimize any downtime. Our team of skilled technicians is trained to diagnose and repair any issues that may arise, ensuring that your machine is always operating at peak performance.

At Picosecond Laser Machine, we are committed to providing our customers with exceptional support and service. Contact us today to learn more about how we can help you get the most out of your investment.

Packing and Shipping:

Product Packaging:

1 Picosecond Laser Machine

1 Power cord

1 Foot pedal

1 Laser protective eyewear

1 User manual

Shipping:

Shipping method: UPS Ground

Estimated delivery time: 3-5 business days

Shipping cost: Free

Tracking number will be provided within 24 hours of shipment

FAQ:

Q: What is the brand name of the Picosecond Laser Machine?

A: The brand name of the Picosecond Laser Machine is LZR.

Q: What is the model number of the Picosecond Laser Machine?

A: The model number of the Picosecond Laser Machine is P6.

Q: Where is the Picosecond Laser Machine manufactured?

A: The Picosecond Laser Machine is manufactured in Guangzhou, China.

Q: What is the wavelength range of the Picosecond Laser Machine?

A: The Picosecond Laser Machine has a wavelength range of 1064nm, 532nm, and 755nm.

Q: What kind of skin concerns can the Picosecond Laser Machine treat?

A: The Picosecond Laser Machine is effective in treating various skin concerns such as pigmentation, tattoo removal, acne scars, and fine lines.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou