




Safety Hair Removal Laser Machine , Permanent Hair Removal Machine For Female

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: TM-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Hair Removal Laser Machine
- Machine Power: 1200W
- Cooling Fan: 2pcs Fans
- Energy Density: Up To 120 J/cm²
- Machine App: Permanent Laser Hair Removal Machines
- Laser Type: Other
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Coolingtype: TEC+Sapphire Contact Cooling
- Style: Stationary, Vertical
- Highlight: **CE Hair Removal Laser Machine,
Safety Hair Removal Laser Machine,
Permanent Hair Removal Machine For Female**



More Images



Product Description

Say Hello To Smooth And Hair-Free Skin With Stationary Laser Hair Elimination Machine

Product Description:

The Hair Removal Laser Machine has a pulse duration of 10-400 Ms, which means that it can deliver precise laser energy to your skin for a short period of time. This ensures that the laser targets the hair follicles without damaging the surrounding skin. Additionally, the machine has an energy density of up to 120 J/cm², which allows it to deliver a high amount of energy to the hair follicles in a short amount of time.

One of the standout features of the Hair Removal Laser Machine is its 15-inch touch screen display. This display is easy to use and allows you to control the settings of the laser with just a few taps of your finger. The touch screen also provides real-time feedback on the energy levels and pulse duration of the laser, so you can adjust the settings as needed to get the best results.

Safety is a top priority with the Hair Removal Laser Machine. The device is equipped with skin contact cooling technology, which cools the skin before and after each laser pulse to minimize discomfort and reduce the risk of skin damage. Additionally, the machine has an automatic power cut-off feature that stops the laser from firing if it is not in contact with the skin.

The Hair Removal Laser Machine also comes with a companion app that allows you to track your hair removal progress and get personalized recommendations for your skin type and hair color. The app is easy to use and provides a wealth of information that can help you get the best possible results from your laser hair removal treatments.

Overall, the Hair Removal Laser Machine is a highly effective and easy-to-use professional hair removal device that can deliver salon-quality results in the comfort of your own home. With its advanced features and safety measures, this machine is a great investment for anyone looking to achieve smooth, hair-free skin.

Features:

Product Name: Hair Removal Laser Machine

Safety Features: Skin Contact Cooling, Automatic Power Cut-off

Laser Type: Other

Style: Stationary, Vertical


Cooling Type: TEC+Sapphire Contact Cooling

Display: 15-inch Touch Screen

This Hair Removal Laser Machine is the perfect Unwanted Hair Remover Device for those looking to achieve smooth skin with laser technology.







755nm for white skin(fine, golden hair)

808nm for yellow/neutral skin

1064nm for black(black hair)

Before and after hair removal



Ice Hair Removal



Selectable
Spot Size



Diode Laser



Fast Light



Adjustable
Frequency



Adjustable Energy



Technical Parameters:

Product Name:	Unwanted Hair Remover Device
Machine App:	Permanent Laser Hair Removal Machines
Energy Density:	Up To 120 J/cm2
Machine Power:	1200W
Style:	Stationary, Vertical
Safety Features:	Skin Contact Cooling, Automatic Power Cut-off
Laser Type:	Other
Power:	600W
Cooling Fan:	2pcs Fans
Skin Types:	I-VI
Display:	15-inch Touch Screen
Product Description:	Advanced Hair Removal System, Permanent Hair Removal Device

Applications:

The LZR Hair Removal Laser MachineTM-1 is a stationary and vertical laser hair elimination machine that is perfect for use in salons, spas, and other professional settings. The machine is designed to be easy to use, with a 15-inch touch screen display that allows users to adjust the settings and monitor the progress of the treatment.

The LZR Hair Removal Laser MachineTM-1 is suitable for all skin types, ranging from I to VI. It uses a high-energy-density laser that can deliver up to 120 J/cm2 of energy, making it one of the most powerful hair removal machines on the market. This laser hair elimination machine is perfect for people who want to remove hair quickly and effectively, without having to endure the pain and discomfort associated with other hair removal methods.

There are many different occasions and scenarios where the LZR Hair Removal Laser MachineTM-1 can be used. Salons and spas can use this hair depilator machine to offer their clients a high-quality, professional hair removal service. The machine is also ideal for use in medical clinics and hospitals, where it can be used to remove hair prior to surgery.

Individuals who want to remove hair from their own homes can also benefit from the LZR Hair Removal Laser MachineTM-1. This laser hair elimination machine is perfect for people who want to remove hair quickly and easily, without having to go to a salon or spa. The machine is easy to use and can be adjusted to suit the user's individual needs.

In conclusion, the LZR Hair Removal Laser MachineTM-1 is a high-quality laser hair elimination machine that is suitable for use in a wide range of scenarios and occasions. Whether you are a professional looking to offer hair removal services or an individual looking for an effective hair removal solution, this machine is sure to meet your needs.

Customization:

Customize your Permanent Hair Removal Device with the following options:

Brand Name: LZR

Model Number: TM-1

Place of Origin: China Guangzhou

Skin Types: I-VI

Power: 600W

Pulse Duration: 10-400 Ms

Coolingtype: TEC+Sapphire Contact Cooling

Safety Features: Skin Contact Cooling, Automatic Power Cut-off

Upgrade your Hair Depilator Machine with the latest Laser Hair Elimination Machine technology from LZR. Our customizable options ensure the perfect fit for your specific needs and preferences.

Support and Services:

The Hair Removal Laser Machine product comes with comprehensive technical support and services to ensure optimal performance and user satisfaction. Our team of experienced technicians is available to assist with any technical issues and provide guidance on proper use and maintenance of the machine. In addition, we offer regular maintenance and calibration services to keep the machine in top condition. Our product also comes with a warranty to provide peace of mind and protection against any defects or malfunctions. For any questions or concerns, please don't hesitate to contact our customer support team.

FAQ:

Q: What is the brand name of the hair removal laser machine?

A: The brand name of the hair removal laser machine is LZR.

Q: What is the model number of the hair removal laser machine?

A: The model number of the hair removal laser machine is TM-1.

Q: Where is the hair removal laser machine manufactured?

A: The hair removal laser machine is manufactured in China Guangzhou.

Q: Is the hair removal laser machine safe to use?

A: Yes, the hair removal laser machine is safe to use. It has been tested and approved by the relevant authorities and meets all safety standards.

Q: How long does it take to see results from using the hair removal laser machine?

A: Results from using the hair removal laser machine may vary depending on the individual. However, most people notice a reduction in hair growth after just a few sessions.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou