



TEC Sapphire Contact Cooling Permanent Hair Removal Equipment 600W

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: TM-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month




Product Specification

- Name: Permanent Hair Removal Equipment
- Coolingtype: TEC+Sapphire Contact Cooling
- Machine App: Permanent Laser Hair Removal Machines
- Power: 600W
- Cooling Fan: 2pcs Fans
- Display: 15-inch Touch Screen
- Pulse Duration: 10-400 Ms
- Energy Density: Up To 120 J/cm2
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Highlight: Permanent Hair Removal Equipment 600W, CE Permanent Hair Removal Equipment, Permanent Hair Removal System 600W



More Images





Product Description

TEC Sapphire Contact Cooling Hair Removal Laser Machine Up To 120 J/Cm2 For I-VI Skin Types

Product Description:

The Smooth Skin Laser Machine uses TEC+Sapphire Contact Cooling to ensure that the skin is always at a safe temperature during use. This cooling system ensures that the hair removal process is painless and comfortable. The machine also has 2pcs fans that further enhance the cooling system, ensuring that the skin is always protected.

The safety features of the Smooth Skin Laser Machine are top-notch. The Skin Contact Cooling feature ensures that the skin is protected during hair removal. The Automatic Power Cut-off feature ensures that the machine automatically shuts down if it senses any abnormal activity or overheating. These safety features ensure that the user is protected from any harm that may arise during use.

The Smooth Skin Laser Machine uses a laser type that is different from traditional laser hair removal machines. This laser type is safer and more efficient, ensuring that the hair removal process is more effective and lasts longer. The machine is designed to cater to all skin types, ensuring that users with different skin types can all enjoy the benefits of the Smooth Skin Laser Machine.

The Smooth Skin Laser Machine is the perfect solution for people looking for a painless and effective hair removal method. The machine is easy to use, and its safety features ensure that users are always protected during use. The machine's ability to cater to all skin types makes it a versatile product that can be used by people of all ages and skin types.

Features:

Product Name: Hair Removal Laser Machine

Coolingtype: TEC+Sapphire Contact Cooling

Machine App: Permanent Laser Hair Removal Machines

Display: 15-inch Touch Screen

Machine Power: 1200W

Skin Types: I-VI


The Hair Removal Laser Machine is an advanced hair removal system that utilizes laser technology to provide efficient and effective hair removal treatment. It features a 15-inch touch screen display, 1200W machine power, and TEC+Sapphire contact cooling for a comfortable and safe experience. Suitable for skin types I-VI, this laser hair treatment device is perfect for permanent hair removal.

Technical Parameters:

Product Attribute	Value
Product	Advanced Hair Removal System
Type	Hair Removal Laser Machine
Laser Type	Other
Cooling Type	TEC+Sapphire Contact Cooling
Cooling Fan	2pcs Fans
Machine Power	1200W
Style	Stationary, Vertical
Power	600W
Machine App	Permanent Laser Hair Removal Machines
Pulse Duration	10-400 Ms
Skin Types	I-VI
Display	15-inch Touch Screen
Device	Permanent Hair Removal Device







755 nm for white skin
(fine, golden hair removal)
808 nm for yellow neutral skin
(brown hair removal)
1064 nm for black skin
(black hair removal)

 Ice Hair Removal	 Selectable Spot Size	 Diode Laser
 Fast Light	 Adjustable Frequency	 Adjustable Energy

Applications:

The LZR TM-1 is suitable for use in a variety of settings, including beauty salons, spas, and medical clinics. It is an ideal choice for individuals who are looking for a professional hair removal device that can provide long-lasting results.

The LZR TM-1 uses advanced laser technology to target hair follicles, effectively eliminating hair growth. It is suitable for use on all skin types, from type I to type VI. The machine's TEC and sapphire contact cooling system ensures a comfortable and pain-free treatment experience for clients.

The LZR TM-1 is an excellent choice for individuals who want to remove unwanted hair from their face, arms, legs, underarms, and bikini line. It is also suitable for treating larger areas of the body, such as the back or chest.

Whether you are a beauty professional looking for a reliable and effective hair removal solution for your clients or an individual looking for a permanent hair removal option, the LZR TM-1 laser hair elimination machine is an excellent choice. With its advanced technology, easy-to-use interface, and professional-grade performance, this laser hair removal machine is sure to satisfy your hair removal needs.

Customization:

Support and Services:

Our Hair Removal Laser Machine product comes with comprehensive technical support and services to ensure your satisfaction with your purchase. Our team of experts is available to assist you with any questions or issues you may encounter with the product. We also offer training and resources to help you properly operate and maintain the machine. Additionally, we provide repair services and replacement parts if needed. Rest assured, we are committed to providing you with the best possible experience with your Hair Removal Laser Machine.

FAQ:

Q: What is the brand of this hair removal laser machine?

A: The brand of this hair removal laser machine is LZR.

Q: What is the model number of this hair removal laser machine?

A: The model number of this hair removal laser machine is TM-1.

Q: Where is this hair removal laser machine manufactured?

A: This hair removal laser machine is manufactured in China Guangzhou.

Q: Is this hair removal laser machine safe to use?

A: Yes, this hair removal laser machine is safe to use as long as you follow the instructions provided in the user manual.

Q: Can this hair removal laser machine be used on all skin types?

A: This hair removal laser machine is suitable for most skin types, but it may not be effective on very light or white hair.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou