



Electric Powered Skin Care Light Therapy Device 110V 220V for Salon

Our Product Introduction

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: SPA-3
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Skin Care Light Therapy Device
- Condition: New
- Led Chip: 5w LED
- Supplier Type: Top 10 Direct Manufacturer From China
- Relative Humidity: $\leq 85\%$
- Battery Capacity: 1600mAh
- Red Light Wavelength: 660 Nm
- Timing Range: Adjustable For 15-60minutes
- Power Supply: Electric-Powered
- Highlight: Skin Care Light Therapy Device 110V,
Electric Powered Skin Care Light Therapy Device
, Skin Light Therapy Device 220V



Product Description

Top 10 Direct From HineSkin Light Therapy Mac 110V 220V

Product Description:

This machine is brand new and comes with a powerful electric-powered supply that ensures optimal performance. The laser output power is ≤ 5 MW, which is safe and effective for use on all skin types. The LED Phototherapy Machine is also equipped with a voltage of 110V and 220V, making it suitable for use in different locations.

The machine is easy to use and comes with a user-friendly interface that allows you to adjust the settings according to your needs. It also comes with a compact design that makes it easy to store and transport. This LED Phototherapy Machine is perfect for use in beauty salons, spas, and clinics.

The LED Phototherapy Machine uses different wavelengths of light to penetrate the skin and stimulate cells that are responsible for producing collagen and elastin. This process helps to plump up the skin, reduce the appearance of fine lines and wrinkles, and improve overall skin health. The machine also helps to reduce inflammation and redness, making it an effective treatment for acne and rosacea. The LED Phototherapy Machine also uses blue light to kill the bacteria that cause acne, making it an effective treatment for acne-prone skin. The machine also uses red light to stimulate the production of collagen and elastin, which helps to reduce the appearance of fine lines and wrinkles.

The LED Phototherapy Machine is safe and pain-free, and it doesn't require any downtime. The machine is also suitable for use on all skin types and can be used on the face, neck, and chest. The machine is also equipped with a feature that allows you to adjust the intensity of the light, making it suitable for use on sensitive skin.

In conclusion, the LED Phototherapy Machine is an excellent investment for anyone looking to improve the overall health and appearance of their skin. This machine is easy to use, safe, and effective, and it comes with a compact design that makes it easy to store and transport. Whether you're looking to reduce the appearance of fine lines and wrinkles or improve the overall texture and tone of your skin, the LED Phototherapy Machine is the perfect solution.

Features:

Product Name: LED Light Therapy Machine

Lamp Power: 192W

Actual Wattage: 300W

Yellow Wavelength: 589 Nm

Condition: New

Power Supply: Electric-Powered

Features:

LED Facial Light Therapy System

LED Therapy Equipment

Photon Light Therapy System

Technical Parameters:

Lamp Power	192W
Red Light Wavelength	660 Nm
Laser Output Power	≤ 5 Mw
Supplier Type	Top 10 Direct Manufacturer From China
Relative Humidity	$\leq 85\%$
Led Quantity	200pcs
Power Supply	Electric-Powered
Timing Range	Adjustable For 15-60minutes
Condition	New
Voltage	110V 220V

This is the technical parameter table of the ShineSkin Light Therapy Machine, also known as the Eczema Treatment LED Machine. It is a professional light therapy machine that features a lamp power of 192W and a red light wavelength of 660 Nm. The laser output power is ≤ 5 Mw and the supplier type is a top 10 direct manufacturer from China. The device can operate in a relative humidity of $\leq 85\%$ and has 200pcs of LED quantity. The power supply is electric-powered and the timing range is adjustable for 15-60 minutes. The device is new and can operate on both 110V and 220V voltage.

Applications:

Salons: The LZR SPA-3 LED Light Therapy Machine is an excellent addition to any salon that offers skin treatments. It can be used to treat acne, wrinkles, and other skin conditions, making it a popular choice among customers.

Home Use: This machine is also suitable for home use. Its lightweight and portable design make it easy to use at home and store away when not in use. It is perfect for those who want professional-grade skin treatments without having to visit a salon.

Spas: Spas can benefit from the LZR SPA-3 LED Light Therapy Machine as well. It can be used to provide relaxing and rejuvenating skin treatments to their customers. The machine's yellow wavelength of 589 Nm is particularly effective in reducing inflammation and promoting circulation.

Professional Use: Dermatologists and other skin care professionals can also use this machine in their practice. Its high lamp power of 192W ensures that it can handle heavy usage and provide consistent results.

Regardless of the occasion or scenario, the LZR SPA-3 LED Light Therapy Machine is a reliable and effective choice for skin treatments. Its battery capacity of 1600mAh ensures that it can be used for extended periods without needing to be recharged. Additionally, its relative humidity of $\leq 85\%$ means that it is suitable for use in a variety of environments.

Choose the LZR SPA-3 LED Light Therapy Machine for your skin care needs and experience the benefits of a professional-grade device in the comfort of your own home or practice.

Customization:





Support and Services:

The LED Light Therapy Machine is a cutting-edge device designed for multiple skin treatments, including anti-aging, acne reduction, and skin rejuvenation. To ensure optimal performance and longevity of the device, we offer comprehensive product technical support and services, including:

- Troubleshooting guidance for any technical issues that may arise
- Maintenance and repair services to keep your device functioning at its best
- Access to product manuals and instructional materials
- Expert advice on how to use the device for maximum effectiveness and safety

Our team of qualified technicians and customer support specialists are dedicated to providing you with the assistance you need to make the most of your LED Light Therapy Machine. Contact us today for more information or to request support.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebain Village, Baiyun District,
Guangzhou