



7D Independent HIFU Beauty Machine With 5mm-25mm Adjustable Treatment Length

Our Product Introduction

For more products please visit us on hifubeauty-machine.com


Basic Information

- Place of Origin: Guangzhou, China
- Brand Name: other
- Certification: CE
- Model Number: hifu-3
- Minimum Order Quantity: 1 Sets
- Price: Discussible
- Packaging Details: Regular Packaging
- Delivery Time: 5-8 working Days
- Payment Terms: L/C, T/T, D/A, D/P, Western Union, MoneyGram
- Supply Ability: 10000 pieces/month



Product Specification

- Input Voltage: AC100V-240V
- Energy Type: HIFU (High Intensity Focused Ultrasound)
- Length: 5.0 - 25mm (1.0mm / Step)
- Lamp Longevity: 300,000 Shots (IPL)
- Specification: Packing , 64 X 81 X 131 Cm
- Energy Output: 0.1J-2.5J (0.1J/step)
- Supplier Info: Manufactueruer
- Spacing: 1.0-10mm(0.5mm/step)
- Highlight: **7D HIFU Beauty Machine ,
Independent HIFU Beauty Machine**



7D Independent Hifu Beauty Machine / Hifu Machine With 5mm-25mm Adjustable Treatment Length



Following the line on one of the cartridges, could be useful tip to design the treatment area.

- 4.5mm probe(lifting): Acting on SAMA fascia layer to Tighten the muscle bottom and achieve firmness
- 3.0mm probe(Tighting): Aiming on the dermis to stimulate continuous regeneration of collagen
- 2.0 mm probe(Wrinkle Erasing): The world's original treatment probe specializing for eyes and deep lines such as forehead, neck and crow's feet
- 1.5mm probe(Toning): Targeting at the epidermis to improve superficial fine lines, skin tone, skin texture, and fine pores




Body treatment site



1.5mm

3.0mm


4.5mm



HIFU MACHINE uses HIGH INTENSITY FOCUS ULTRASOUND TECHNOLOGY, so that the skin tissue creates the thermal, and make high-speed friction cells to stimulate collagen. Such heat effect will not hurt the epidermis, because the treatment is fast and direct access to the treatment site within 0-0.5 seconds, without touching the surrounding tissue and can be passed directly to the superficial tendon membrane system (SMAS).


THE HIFU MACHINE can be done firm skin while pulling the muscle layer, thin face up to the progressive effect. What is SMAS? Lift operation at present by plastic surgeons do operation layer, the SMAS layer, which is (superficial musculoaponeurotic system, referred to as the fascia (fascia)). SMAS layer on the skin depth is about 4.5mm, subcutaneous fat and muscle.

Safety Is Our Highest Priority



PROCEDURE SAFETY

Ultrasound energy equipped multi-depth cartridges of the 7D HIFU only target specific layers of skin without damaging nearby tissue.



CLINICIAN SAFETY

The 7D HIFU is performed by laboriously-trained practitioners that can assure a professional standard of care and clinical expertise for safe procedures.

Reform Your Youth

The 7D HIFU complies with both clinical and international regulatory standards to assure safe and effective procedures for lifting and tightening the face as well as slimming and tightening the contours of the body.



Products real shot





Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,
Guangzhou