



## Vertical HIFU Beauty Machine 0.1J-2.5J Energy 5-25mm Scope

### Our Product Introduction

#### Basic Information

- Place of Origin: Guangzhou, China
- Brand Name: other
- Certification: CE
- Model Number: hifu-3
- Minimum Order Quantity: 1 Sets
- Price: Discussible
- Packaging Details: Regular Packaging
- Delivery Time: 5-8 working Days
- Payment Terms: L/C, T/T,
- Supply Ability: 10000 pieces/month



#### Product Specification

- Body Type: Vertical
- Specification: Packing , 64 X 81 X 131 Cm
- Energy Type: HIFU (High Intensity Focused Ultrasound)
- Lamp Longevity: 300,000 Shots (IPL)
- Supplier Info: Manufactueruer
- Size: Portable
- Scope: 5.0mm-25mm
- Pack Size: 60\*44\*46cm
- Highlight: **HIFU beauty machine with 300 lamp longevity ,  
HIFU skin tightening machine 1.0-10mm  
spacing  
, HIFU facial machine 0.1J-2.5J energy output**

### Product Description:

The HIFU Beauty Machine, also known as the High Intensity Focused Ultrasound Beauty Machine, is an advanced beauty treatment device designed to provide intense ultrasonic treatment for various skin concerns. Manufactured by a reputable supplier, this cutting-edge machine offers a wide scope of treatment depths ranging from 5.0mm to 25mm, making it suitable for addressing different skin layers with precision.

One of the key features of this HIFU Beauty Machine is its energy output flexibility, which ranges from 0.1J to 2.5J with a step increment of 0.1J. This customizable energy output allows for personalized treatment settings based on individual skin needs and sensitivity levels, ensuring optimal results and safety during the procedure.

The body type of this HIFU Beauty Machine is vertical, providing a stable and ergonomic design for easy operation and maneuverability during treatments. Its vertical orientation also helps save space in the treatment room, making it a practical choice for beauty clinics and aesthetic centers with limited space.

With a spacing range of 1.0mm to 10mm and a step increment of 0.5mm, this HIFU Beauty Machine offers precise control over the treatment area and depth, allowing practitioners to target specific skin concerns with accuracy. Whether performing face lifting treatments or addressing other facial rejuvenation needs, the adjustable spacing feature enhances treatment efficacy and customization.

Overall, the HIFU Beauty Machine is a versatile and efficient ultrasonic beauty system that delivers exceptional results for various skin rejuvenation treatments. From face lifting and tightening to wrinkle reduction and skin firming, this advanced device offers a comprehensive solution for clients seeking non-invasive yet effective beauty enhancements.

### Features:

Product Name: HIFU-Beauty Machine

Input Voltage: AC100V-240V

Specification: Packing, 64 X 81 X 131 Cm

Energy Output: 0.1J-2.5J (0.1J/step)

Energy Type: HIFU (High Intensity Focused Ultrasound)

Length: 5.0 - 25mm (1.0mm/Step)

### Technical Parameters:

Scope	5.0mm-25mm
Length	5.0 - 25mm (1.0mm / Step)
Body Type	Vertical
Specification	Packing, 64 X 81 X 131 Cm
Supplier Info	Manufacturer
Spacing	1.0-10mm (0.5mm/step)
Lamp Longevity	300,000 Shots (IPL)
Energy Output	0.1J-2.5J (0.1J/step)
Pack Size	60*44*46cm
Material	ABS

### Applications:

The other HIFU-3 Non-invasive HIFU Skin Tightening Device is a cutting-edge beauty machine designed for a variety of application occasions and scenarios. With its advanced technology and exceptional features, this Non-invasive Ultrasound Beauty Therapy Machine is perfect for professionals in the beauty industry looking to offer top-notch services to their clients.

One of the key product attributes of the HIFU-3 is its versatile working mode, which allows for both simultaneous and independent control. This feature makes it suitable for a wide range of beauty treatments, including face lifting, skin tightening, and overall rejuvenation.

Whether used in a beauty salon, spa, or cosmetic clinic, the HIFU-3 delivers impressive results and ensures client satisfaction. Its Scope of 5.0mm-25mm makes it ideal for targeting different skin depths, providing customized treatments for various skin concerns.

Manufactured in Guangzhou, China, this CE-certified Ultrasonic Beauty System for Face Lifting guarantees high quality and performance. The product's ABS material and compact Pack Size of 60\*44\*46cm make it durable and easy to transport, while the regular Packaging Details ensure safe delivery.

With a Minimum Order Quantity of 1 set and a discussible Price, the HIFU-3 is accessible to a wide range of beauty professionals. The fast Delivery Time of 5-8 working days, flexible Payment Terms (L/C, T/T), and impressive Supply Ability of 10000 pieces per month make it a convenient and reliable choice for businesses.

In conclusion, the HIFU-3 Non-invasive HIFU Skin Tightening Device is a must-have tool for professionals seeking effective and efficient beauty treatments. Its innovative technology, versatile working modes, and exceptional results make it a valuable asset in any beauty setting.



## Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubeaauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou